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21.1 ใบแทรกน้ำยา Creatinine Flex® reagent cartridge (PI-LAB-003)

21.2 SOP For Dimension EXL 200 Clinical Chemistry System (MN-LAB-002)

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REF DF270B

PI-LAB-003/00(01/10/2560)

Dimension® clinical chemistry system

Flex® reagent cartridge

EZCR

Issue Date 2011-03-07

Enzymatic Creatinine

Intended Use: The EZCR method is an *in vitro* diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension® clinical chemistry system. Creatinine measurements are used in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for other urine analytes.

Summary: This creatinine method measures creatinine enzymatically. Enzymatic creatinine methods are reported to be less susceptible to interference than Jaffe methods from non-creatinine substances.¹


Serum and plasma creatinine concentrations, freely filtered by the glomerulus, are inversely related to glomerular filtration rate (GFR).² Therefore, both blood and urine creatinine measurements are used to assess kidney function and in the diagnosis and treatment of renal diseases. Creatinine is also useful in monitoring effectiveness of renal dialysis.

Additionally, calculated creatinine clearance values may be used to establish effective therapeutic dosing levels of pharmaceuticals.³

Test results of the Dimension® EZCR method correlate with those obtained by an IDMS-traceable method.

Principles of Procedure: In a TAPS** buffered coupled enzyme reaction, creatininase hydrolyzes creatinine to creatine, which is hydrolyzed by creatinase to sarcosine. Sarcosine oxidase hydrolyzes sarcosine to glycine, formaldehyde, and peroxide. The peroxide and a chromogen in the presence of peroxidase form a colored end product that is proportional to the amount of creatinine in the sample and is measured using a bichromatic (540, 700 nm) end point technique.

**3-[[tris(hydroxymethyl)methyl]amino]propanesulfonic acid (TAPS)

*HTIB = 2,4,6-triiodo-3-hydroxybenzoic acid

Reagents

Wells*	Form	Ingredient	Concentration*	Source
1 – 2	Liquid	HCl	0.5N	
3 – 4 ^c	Liquid	TAPS buffer	30 mmol/L	Bacterial
		Creatininase	>333 µkat/L	
		Sarcosine oxidase	>33 µkat/L	
		HTIB	5.9 mmol/L	
5 – 6 ^c	Liquid	TAPS buffer	50 mmol/L	Bacterial
		Creatininase	>500 µkat/L	
		Horseradish peroxidase	≥16.7 µkat/L	
		4-aminophenazone	2.0 mmol/L	
		Potassium hexacyanoferrate (II)	163 µmol/L	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Wells 3 – 6 contain detergents and preservatives.

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1 – 6

Specimen Collection and Handling: Recommended specimen types: serum, sodium and/or lithium heparin plasma, and urine.

Grossly hemolyzed samples should not be used with the EZCR method.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.⁴

Follow the instructions provided with your specimen collection device for use and processing.⁵

For serum, complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁶ Separated serum/plasma samples should be stored at 2 – 8 °C and analyzed within 48 hours. For longer storage, samples should be frozen at or below -20 °C.⁸

Urine creatinine measurements require no special patient preparation for specimen collection.

Collect random and timed urine specimens using recommended procedures for collection, transportation, and preservation of urine specimens.⁷

Random or 24 hour urine collections require no preservatives. Specimens previously preserved with 6N HCl or Boric Acid are acceptable.⁷

Urines (random or 24 hour collections) should be stored at 2 – 8 °C and analyzed within 4 days. Freeze for longer storage.⁸

Urine specimens should be free of particulate matter.

The purpose of specimen storage information is to provide guidance to users; however, users may validate their own procedures for storing patient samples.

Procedure

Materials Provided

EZCR Flex® reagent cartridge, Cat. No. DF270B

Materials Required But Not Provided

Dimension® CHEM I Calibrator, Cat. No. DC18B

Quality control materials

Test Steps

Sampling⁹, reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® clinical chemistry system manual.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	6 µL
Reagent 1 Volume	110 µL
Reagent 2 Volume	118 µL
Temperature	37.0 °C ± 0.1 °C
Reaction Time	8.7 minutes ⁸
Wavelength	540 and 700 nm
Type of Measurement	Bichromatic end point
e. calculated as time from test initiation to final result	

Calibration

Measurement Range 0.03 – 20.00 mg/dL [3 – 1768 µmol/L]¹
 Calibration Material Dimension® CHEM I Calibrator Cat. No. DC18B
 Calibration Scheme 3 levels in triplicate
 Units mg/dL [µmol/L]
 (mg/dL x 88.4) = [µmol/L]

Typical Calibration Levels

Level 1: 0.00 mg/dL [0 µmol/L]
 Level 2: 11.00 mg/dL [972 µmol/L]
 Level 3: 22.00 mg/dL [1945 µmol/L]

Calibration Frequency

A new calibration is required

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

Assigned Coefficients
 C₀ -0.210
 C₁ 0.070

f. Système International d'Unités [SI Units] are in brackets.

Quality Control

Follow government regulations or accreditation requirements for quality control frequency.

At least once each day of use, analyze two levels of a Quality Control (QC) material with known creatinine concentrations. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument calculates the concentration of creatinine in mg/dL [µmol/L] using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR) for serum/plasma: 0.03 – 20.00 mg/dL [3 – 1768 µmol/L]

Analytical Measurement Range (AMR) for urine: 0.6 – 400.00 mg/dL [53 – 35360 µmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process.

• **Serum/plasma samples** with results in excess of 20.00 mg/dL [1768 µmol/L] are reported as "Above Assay Range" and should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water to obtain results within the analytical measurement range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The autodilute sample volume is 2 µL (dilution factor = 3) for serum/plasma. Refer to your Dimension® Operator's Guide.

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- Automated Urine Dilution (AUD):** Urine samples are automatically diluted 1:20 with the Dimension® System water (1 part urine sample and 19 parts water) resulting in an analytical measurement range of 0.6 – 400.00 mg/dL [53 – 35360 µmol/L]. Urine samples with results in excess of 400.00 mg/dL [35360 µmol/L] are reported as "Above Assay Range" and should be repeated on dilution.
Manual Dilution: Dilute with Reagent grade water to obtain results within the analytical measurement range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.
Autodilution (AD) is not available for urines that have the AUD feature.
- Serum/plasma** with results less than 0.03 mg/dL [3 µmol/L] should be reported as "less than 0.03 mg/dL [3 µmol/L]". **Urines** with results less than 0.6 mg/dL [53 µmol/L] should be reported as "less than 0.6 mg/dL [53 µmol/L]".

Calculated result applications available in the Dimension® EasyLink™ Informatics System:

Creatinine clearance (mL/min)^{8,9} = (U (mg/dL) x (V (mL/min))/(S (mg/dL)) where

U = urine creatinine in mg/dL

V = urine volume in mL/min

S = serum creatinine in mg/dL

An estimated glomerular filtration rate (eGFR) can be calculated using Dimension® EZCR results and one of the following Modification of Diet in Renal Disease (MDRD) Study equations:

IDMS-traceable MDRD Study Equation eGFR calculation¹⁴ for use with EZCR method:

eGFR (mL/min/1.73 m²) = **175** x (EZCR mg/dL)^{-1.154} x (Age)^{-0.203} x (0.742 if female) x (1.210 if African American) (conventional units)

eGFR (mL/min/1.73 m²) = **175** x (EZCR µmol/L/88.4)^{-1.154} x (Age)^{-0.203} x (0.742 if female) x (1.210 if African American) [SI units]

Use of equations such as those shown above to report estimated GFR should be limited to adults between the ages of 18 and 75 years who are of average body composition and are not experiencing rapidly changing kidney function status.¹⁰

Estimation of GFR using mathematical equations such as the MDRD eGFR is NOT recommended when patients have: extremes of age or body size, severe malnutrition or obesity, skeletal muscle diseases, paraplegia or quadriplegia, limb amputations, vegetarian diets, pregnancy, rapidly changing kidney function, prior to dosing drugs with significant toxicity that are excreted by kidneys, or any other condition in which the patient does not have stable kidney function. Such patients may need measurement of clearance instead of estimation.¹⁰

Clinical practice guidelines with limitations for use of the MDRD Study equation for reporting eGFR can be found on the NCKDE, KDOQI, and IFKF websites.^{10,11,12} Any future refinements in the recommended MDRD Study equation(s) would also be published there.

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in EZCR results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported if indicated. A system malfunction may exist if the following 5-test precision is observed:

EZCR Concentration	SD
1.00 mg/dL [88.4 µmol/L]	>0.03 mg/dL [2.65 µmol/L]
6.00 mg/dL [530.4 µmol/L]	>0.18 mg/dL [15.9 µmol/L]

Interfering Substances

The EZCR was evaluated for interference according to CLSI/NCCLS EP7-A2.¹³ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Hemoglobin (hemolysate) at 400 mg/dL [0.25 mmol/L] increases EZCR results by 15% at a creatinine concentration of 1.00 mg/dL [88.4 µmol/L].

Blood samples from patients with Waldenström's macroglobulinemia may contain monoclonal IgM that produce falsely elevated results.¹⁴

Expected Values

Serum/plasma¹⁵

Males: 0.67 – 1.17 mg/dL [59 – 104 µmol/L]

Females: 0.51 – 0.95 mg/dL [45 – 84 µmol/L]

Random Urine¹⁵

Males: 40 – 278 mg/dL [3536 – 24575 µmol/L]

Females: 29 – 226 mg/dL [2564 – 19978 µmol/L]

24 Hour Urine Excretion¹⁵

Males: 0.87 – 2.41 g/day [8 – 21 mmol/day]

Females: 0.67 – 1.59 g/day [6 – 14 mmol/day]

Each laboratory should establish its own expected values for EZCR as performed on the Dimension® clinical chemistry system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® clinical chemistry system.

Material	Precision ^{17,a}		Standard Deviation (% CV)	
	Mean mg/dL [µmol/L]	Repeatability	Within-Lab	
			Repeatability	Within-Lab
Serum				
Serum pool 1	1.04 [92]	0.01 [0.88] (1.1)	0.02 [1.77] (2.0)	
Serum pool 2	1.55 [137]	0.01 [0.88] (0.7)	0.02 [1.77] (1.5)	
BioRad® Liquichek™				
Level 1	0.79 [70]	0.01 [0.88] (1.3)	0.03 [2.65] (3.4)	
Level 2	5.80 [512]	0.03 [2.65] (0.5)	0.05 [4.42] (0.9)	
Urine				
Urine pool 1	116.77 [10322]	1.10 [97] (0.9)	1.73 [153] (1.5)	
Urine pool 2	293.38 [25935]	2.82 [249] (1.0)	4.02 [355] (1.4)	
MAS® Urine Chemistry Liquid				
Level 1	66.94 [5917]	0.94 [83.1] (1.4)	1.08 [95.5] (1.6)	
Level 2	139.42 [12325]	1.61 [142.3] (1.2)	1.80 [159.1] (1.3)	

g. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Liquichek™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

MAS® is a registered trademark of Microgenics Corporation, Fremont, CA 94538, USA

Method Comparison¹⁸

Regression Statistics^a

Comparative Method	Slope	Intercept mg/dL [µmol/L]	Correlation Coefficient	n
Serum/Plasma				
Roche/Hitachi CREA-plus	1.03	-0.03 [-2.7]	1.000	123 ⁱ
Urine				
Roche/Hitachi CREA-plus	0.97	+1.63 [144]	0.998	63 ^j

h. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was orthogonal regression fit.

i. The range of serum and plasma creatinine values in the correlation study was 0.64 – 18.95 mg/dL [56.6 – 1675 µmol/L].

j. The range of urine creatinine values in the correlation study was 11.9 – 400.00 mg/dL [1052 – 35360 µmol/L].

Roche/Hitachi CREA-plus is a product of Roche Diagnostics Division, GmbH, D68298 Mannheim.

Specificity

Hemolysis, Icterus, Lipemia (HIL) Interference

The EZCR was evaluated for interference according to CLSI/NCCLS EP7-A2.¹³ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Creatinine mg/dL [µmol/L]	Bias* %
Hemoglobin (hemolysate)	300 mg/dL [0.19 mmol/L]	1.00 mg/dL [88.4 µmol/L]	<10
	400 mg/dL [0.25 mmol/L]	1.00 mg/dL [88.4 µmol/L]	+15
Bilirubin (unconjugated)	23 mg/dL [393 µmol/L]	1.00 mg/dL [88.4 µmol/L]	<10
	20 mg/dL [342 µmol/L]	1.00 mg/dL [88.4 µmol/L]	<10
Lipemia (Intralipid®)	3000 mg/dL [33.9 mmol/L]	1.22 mg/dL [108 µmol/L]	<10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* **Analyte results should not be corrected based on this bias.**

Non Interfering Substances – serum/plasma

The following substances do not interfere with the EZCR method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at EZCR method concentrations of 1.00 mg/dL [88.4 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20.0 mg/dL	1324 µmol/L
Amikacin	8.0 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	6.0 mg/dL	342 µmol/L
Benzyl alcohol (in heparin)	1%	1%
Caffeine	6.0 mg/dL	308 µmol/L
Carbamazepine	3.0 mg/dL	127 µmol/L
Cefazolin	1.5 mg/dL	133 µmol/L
Cephalexin	1.5 mg/dL	133 µmol/L
Cephalexidine	1.5 mg/dL	133 µmol/L
Cephalexin	1.5 mg/dL	133 µmol/L
Cephadrine	1.5 mg/dL	133 µmol/L
Chloramphenicol	5.0 mg/dL	155 µmol/L
Chloridiazepoxide	1.0 mg/dL	33.3 µmol/L
Chlorpromazine	0.20 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2.0 mg/dL	79.2 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	18 µmol/L
Digoxin	5.0 ng/mL	7.8 nmol/L
Dipyron (metamizol)	5.0 mg/dL	150 µmol/L
D/L-proline	11.5 mg/dL	1017 µmol/L
Dobutamine	0.5 mg/dL	44.2 µmol/L
Dopamine	0.1 mg/dL	6.5 µmol/L
Erythromycin	6.0 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25.0 mg/dL	1770 µmol/L
Furosemide	6.0 mg/dL	181 µmol/L
Gentamicin	1.0 mg/dL	21 µmol/L
Heparin	3.0 U/mL	3000 U/L
Ibuprofen	50.0 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Immunoglobulin M (IgM)	10 g/L	10 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
N-acetyl-cysteine	45 mg/dL	276 µmol/L
N-ethylglycine (lidocaine metabolite)	5.0 mg/dL	442 µmol/L
Nicotine	0.10 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8.0 mg/dL	354 µmol/L
Phenobarbital	10.0 mg/dL	431 µmol/L
Phenytoin	5.0 mg/dL	198 µmol/L
Primidone	4.0 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4.0 mg/dL	222 µmol/L
Urea	500 mg/dL	83 mmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Vancomycin	10.25 mg/dL	71 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

**Non Interfering Substances – urine**

The following substances do not interfere with the EZCR method when present in urine at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at EZCR method concentrations of 45.00 mg/dL [3978 µmol/L].

Substance	Test Concentration	SI Units
5% NaOH	0.02%	0.02%
6N HCl	0.01%	0.01%
Acetone	1.0 g/dL	1.72 mol/L
Ascorbic Acid	1.5 g/dL	45.4 µmol/L
Bilirubin, conjugated	2.0 mg/dL	34.2 µmol/L
Bilirubin, unconjugated	2.0 mg/dL	34.2 µmol/L
Boric acid	1% w/v	162 mmol/L
Cefazolin	1.5 mg/dL	133 µmol/L
Cephalexin	1.5 mg/dL	133 µmol/L
Cephalexidine	1.5 mg/dL	133 µmol/L
Cephalexin	1.5 mg/dL	133 µmol/L
Cephadrine	1.5 mg/dL	133 µmol/L
D/L-proline	11.5 mg/dL	1017 µmol/L
Ethanol	1.0 g/dL	217 mmol/L
Gamma globulin	0.5 g/dL	5.0 g/L
Glucose	2.0 g/dL	0.11 mol/L
Hemoglobin	115 mg/dL	0.07 mmol/L
Human serum albumin	0.5 g/dL	5.0 g/L
N-ethylglycine	10 mg/dL	884 µmol/L
Oxalic acid	0.1 g/dL	110 mmol/L
pH Effects	1.2 to 12.4	1.2 to 12.4
Riboflavin	7.5 mg/dL	2.0 mmol/L
Sodium azide	1% w/v	154 mmol/L
Sodium chloride	6.0 g/dL	10.3 mol/L
Sodium fluoride	1% w/v	238 mmol/L
Urea	6.0 g/dL	10 mol/L

Cross-reactivity:

Creatine and pyruvate were evaluated for cross-reactivity in the presence of 1.0 mg/dL [88.4 µmol/L] of creatinine. The percent cross-reactivity was calculated as follows:

$$\% \text{ cross-reactivity} = \frac{\text{measured analyte} - \text{control analyte (x units [SI])}}{\text{metabolite added}} \times 100$$

Substance	Unit	SI	% cross-reactivity
Creatine	25 mg/dL	1907 µmol/L	5
Pyruvate	6 mg/dL	684 µmol/L	3

Limit of detection¹⁵: 0.03 mg/dL [3 µmol/L]

The limit of detection represents the lowest concentration of creatinine that can be detected with at least 95% probability. CLSI/NCCLS EP17-A, Protocols for Determination of Limits of Detection and Limits of Quantitation, was followed using 688 replicate determinations at two concentration levels to determine Limit of Blank and Limit of Detection.


Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex® and EasyLink™ are trademarks of Siemens Healthcare Diagnostics.

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All rights reserved.