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Material	Precision ^{12,e,f}		
	Mean mg/dL [μmol/L]	Standard Deviation (% CV)	
		Repeatability SD mg/dL [μmol/L]	Within-Lab SD mg/dL [μmol/L]
Serum pool 1	12.2 [212]	0.04 [0.7] (0.4)	0.11 [2] (0.9)
Serum pool 2	21.5 [368]	0.08 [1] (0.4)	0.16 [3] (0.8)
BioRad Lyphocheck® Control Level 1	0.84 [14]	0.01 [0.2] (1.7)	0.02 [0.3] (2.7)
BioRad Lyphocheck® Control Level 2	4.3 [74]	0.02 [0.3] (0.4)	0.04 [0.7] (1.0)
MAS® Liquid Bilirubin QC Level 3	18.5 [316]	0.08 [1] (0.4)	0.2 [3] (1.0)

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Lyphocheck® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618.

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Method Comparison

Comparative Method	Regression Statistics ^g			Correlation Coefficient	n ^{h,i}
	Slope	Intercept mg/dL	SI Units		
Dimension® TBIL	0.937	-0.058	[-1]	1.000	157
Vitros™ TBIL	1.024	-0.337	[-6]	0.996	169

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REF DF167

PI-LAB-016/00(01/10/2560)

Dimension® clinical chemistry system

Flex® reagent cartridge

TBI

See shaded sections: Updated information from 2007-01 version.

Issue Date 2008-04-07

Total Bilirubin

Intended Use: The TBI method for the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended to quantitatively measure total bilirubin in human serum and plasma. Measurements of total bilirubin are used in the diagnosis and treatment of liver, hemolytic, hematological, and metabolic disorders, including hepatitis and gallbladder disease.

Summary: There are at least four distinct bilirubin fractions that make up total bilirubin in serum. The direct reacting fractions are mono- and diconjugated bilirubin (δ and γ-bilirubin) and the delta fraction (δ-bilirubin), which is tightly bound to albumin. Unconjugated bilirubin (α-bilirubin) is water-insoluble and reacts only after addition of an accelerator such as caffeine.¹ The TBI method is a modification of the Doumas reference method,² which is a modification of the diazo method described by Jendrasik and Grof in 1938.³

Principles of Procedure: Diazotized sulfanilic acid is formed by combining sodium nitrite and sulfanilic acid at low pH. Bilirubin (unconjugated) in the sample is solubilized by dilution in a mixture of caffeine/benzoate/acetate/EDTA. Upon addition of the diazotized sulfanilic acid, the solubilized bilirubin including conjugated bilirubins (mono and diglucuronides) and the delta form³ (biliprotein-bilirubin covalently bound to albumin) is converted to diazo-bilirubin, a red chromophore representing the total bilirubin which absorbs at 540 nm and is measured using a bichromatic (540, 700 nm) endpoint technique. A sample blank correction is used.

Solubilized bilirubin + Diazotized sulfanilic acid → Red chromophore (absorbs at 540 nm)

Reagents

Wells*	Form	Ingredient	Concentration*
1, 4 – 6	Liquid	Acetate Buffer	
		Caffeine	168 mM
		Sodium Benzoate	337 mM
		Disodium EDTA	2.57 mM
2	Liquid	Sulfanilic acid	25.89 mM
		Hydrochloric acid	132 mM
		Sodium Nitrite	72.5 mM
3	Liquid		

- a. Wells are numbered consecutively from the wide end of the cartridge.
b. Nominal value per well in a cartridge.

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 5 days for wells 1, 4 – 6
3 days for well 3
15 days for well 2 (diazotized sulfanilic acid formed by the automatic addition of sodium nitrite from well 3)

Specimen Collection and Handling: Serum, lithium heparin plasma, and EDTA plasma can be collected by normal procedures.^{3, 8}

Follow the instructions provided with your specimen collection device for use and processing.⁷

Serum and plasma specimens should be separated from cells within 2 hours after venipuncture.⁸

Specimens should be free of particulate matter. To prevent the appearance of fibrin in serum samples, complete clot formation should take place before centrifugation.⁸ Clotting time may be increased due to thrombolytic or anticoagulant therapy.

Bilirubin is photosensitive. Care should be taken to protect sample from both daylight and fluorescent light to avoid photodegradation. Separated specimens are stable for 8 hours at room temperature, 7 days at 2 – 8 °C or 6 months frozen at -20 °C or colder. Protection from light is required when specimens are stored for more than 8 hours.⁹

Procedure

Materials Provided

TBI Flex® reagent cartridge, Cat. No. DF167

Materials Required But Not Provided

TBI/DBI Calibrator Cat. No. DC167

Quality Control Materials

Purified Water Diluent (Cat. No. 710615901) or Reagent grade water

Test Steps

Sampling,⁴ reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

- c. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

	Cuvette
Sample Volume	10 µL
Reagent 1 Volume	250 µL
Reagent 2 Volume	47 µL
Temperature	37 °C
Wavelength	540, 700 nm
Type of Measurement	Bichromatic endpoint

Calibration

Assay Range	0.1 – 25.0 mg/dL [2 – 428 µmol/L] ¹
Calibration Material	TBI/DBI Calibrator, Cat. No. DC167
Calibration Scheme	3 levels, n = 3
Units	mg/dL [µmol/L] (mg/dL x 17.1) = [µmol/L]
Typical Calibration Levels	0.00, 10.00, 25.00 mg/dL 0, 171, 428 µmol/L
Calibration Frequency	Every 90 days for any one lot
A new calibration is required	<ul style="list-style-type: none"> For each new lot of Flex® reagent cartridges After major maintenance or service, if indicated by quality control results As indicated in laboratory quality control procedures When required by government regulations
Assigned Coefficients	C ₀ 0.00 C ₁ 0.078

Note: Level 1 calibrator for TBI is not included in the TBI/DBI calibrator carton. Purified Water Diluent (Cat. No. 710615901) or Reagent grade water should be used as the Level 1 calibrator for the TBI method.
d. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known total bilirubin concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of total bilirubin in mg/dL [µmol/L] using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0.1 – 25.0 mg/dL [2 – 428 µmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 25.0 mg/dL [428 µmol/L] should be repeated on dilution.

Manual Dilution: Dilute with Reagent grade water 1:4 (1 part of sample with 3 parts of Reagent grade water) to obtain results within the assay range. Enter dilution factor (4). Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The recommended autodilute sample volume is 5 µL for serum and plasma. Refer to your Dimension® Operator's Guide.

Results less than 0.1 mg/dL [2 µmol/L] should be reported as "less than 0.1 mg/dL [2 µmol/L]" .

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed:

TBI Concentration	SD
1.1 mg/dL [19 µmol/L]	> 0.03 mg/dL [0.5 µmol/L]
18.8 mg/dL [322 µmol/L]	> 0.56 mg/dL [10 µmol/L]

Interfering Substances

The TBI method was evaluated for interference according to CLSI/NCCLS EP7-A.¹⁰ Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in mg/dL [µmol/L]. Bias exceeding 10% is considered interference.

Interferent	Concentration SI Units	Total Bilirubin	Bias	Bias (%)
Levodopa	300 µg/mL	1.1 mg/dL	+2.6 mg/dL	+236
	[1.52 mmol/L]	[19 µmol/L]	[44 µmol/L]	
Lipemia (Intralipid®)	600 mg/dL	1.1 mg/dL	+0.2 mg/dL	+18
	[6.78 mmol/L]	[19 µmol/L]	[3 µmol/L]	
Phenazopyridine	80 µg/mL	1.1 mg/dL	+0.8 mg/dL	+42
	[320 µmol/L]	[19 µmol/L]	[14 µmol/L]	

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

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Expected Values: 0.2 – 1.0 mg/dL [3 – 17 µmol/L]¹¹

The reference interval was validated in a confirmatory study using 30 serum samples.

Each laboratory should establish its own reference interval for total bilirubin as performed on the Dimension® system.

Specific Performance Characteristics

All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

Material	Precision ^{12,13}		
	Standard Deviation (% CV)		
	Mean mg/dL [µmol/L]	Repeatability SD mg/dL [µmol/L]	Within-Lab SD mg/dL [µmol/L]
Serum pool 1	12.2 [212]	0.04 [0.7] (0.4)	0.11 [2] (0.9)
Serum pool 2	21.5 [368]	0.08 [1] (0.4)	0.16 [3] (0.8)
BioRad Lyphochek® Control Level 1	0.84 [14]	0.01 [0.2] (1.7)	0.02 [0.3] (2.7)
BioRad Lyphochek® Control Level 2	4.3 [74]	0.02 [0.3] (0.4)	0.04 [0.7] (1.0)
MAS® Liquid Bilirubin QC Level 3	18.5 [316]	0.08 [1] (0.4)	0.2 [3] (1.0)

e. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Quantitative Measurement Methods (EP5-A2, 2004).

f. Specimens at each level were analyzed in duplicate, twice a day, for 20 days. The repeatability (within-run) and within-lab (total) standard deviations were calculated by analysis of variance method.

Lyphochek® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618.

MAS® is a registered trademark of Medical Analysis Systems, Inc., Camarillo, CA 93012.

Method Comparison

Comparative Method	Regression Statistics ¹⁴				
	Slope	Intercept mg/dL	SI Units	Correlation Coefficient	n ¹⁵
Dimension® TBL	0.937	-0.058	[-1]	1.000	157
Vitros™ TBL	1.024	-0.337	[-6]	0.996	169

g. Model equation for regression statistics is: results for Dimension® system = [slope x (comparative method results)] + intercept.

h. The range of TBL values in the Dimension® correlation study was: 0.2 – 22.4 mg/dL [3 – 383 µmol/L].

i. The range of TBL values in the Vitros™ TBL correlation study was: 0.2 – 24.9 mg/dL [3 – 426 µmol/L].

Vitros™ is a trademark of Ortho Clinical Diagnostics, Inc., Rochester, NY 14626.

Recommended Sample Types

Matched human serum and plasma samples were analyzed for total bilirubin. As shown in the table below, no clinically significant difference was observed between different sample types using ordinary least squares analysis to fit the regression line.

Sample Type	n	Slope	Intercept (b)	Correlation Coefficient (r)
Lithium heparin plasma versus serum	25	1.01	-0.11	0.998
EDTA plasma versus serum	25	1.03	-0.30	0.997
EDTA plasma versus lithium heparin plasma	25	0.98	0.13	0.998

Hemolysis and Lipemia Interference Information

The TBL method was evaluated for interference according to CLSI/NCCLS EP7-A. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Test Concentration	TBL Concentration	Bias ¹
	SI Units	SI Units	%
Hemoglobin ¹⁶	1000 mg/dL	1.1 mg/dL	<10
	[0.62 mmol/L] (monomer)	[19 µmol/L]	
	1000 mg/dL	24.8 mg/dL	<10
Lipemia (Intralipid®) ¹⁷	[0.62 mmol/L] (monomer)	[424 µmol/L]	
	200 mg/dL	1.1 mg/dL	<10
	[2.26 mmol/L]	[19 µmol/L]	
	600 mg/dL	1.1 mg/dL	+18
	[6.78 mmol/L]	[19 µmol/L]	

j. Analyte results should not be corrected based on this bias.

k. Samples with hemolysis greater than 1000 mg/dL [0.62 mmol/L] of hemoglobin will be flagged with a "Hemoglobin" error message. Refer to your Dimension® Operator's Guide.

l. Samples with Intralipid® concentrations greater than 600 mg/dL [6.78 mmol/L] may be flagged with an "Abnormal Reaction" error message. Refer to your Dimension® Operator's Guide.

Non-Interfering Substances

The following substances do not interfere with the TBL method when present in serum in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at a total bilirubin concentration of 1.0 ± 0.4 mg/dL [17 ± 7 µmol/L] and 25.0 ± 2.0 mg/dL [428 ± 34 µmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1328 µmol/L
Amikacin	15 mg/dL	256 µmol/L
4-aminosalicylic acid	80 µg/dL	5 µmol/L
Amphotericin B	4 µg/dL	4 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Biliverdin	4 mg/dL	0.061 mmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Carbenicillin	3 mg/dL	7.1 mmol/L
Cefotiam	840 µg/mL	1.4 mmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	1 mg/dL	31.4 µmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40 ¹⁸	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Digoxin	5 ng/mL	6.15 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5 g/dL	50 g/L
Levodopa ¹⁹	300 µg/mL	1.52 mmol/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Methotrexate	50 µg/mL	110 µmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nitrofurantoin	20 µg/mL	0.08 mmol/L
Oxyltetraacycline	50 mg/dL	0.10 mmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	10 mg/dL	443 µmol/L
Phenazopyridine ²⁰	80 µg/mL	320 µmol/L
Phenobarbital	15 mg/dL	647 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Pindolol	5 mg/L	0.02 mmol/L
Piroxicam	10 µg/mL	30 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	1 mg/dL	25 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Reserpine	60 mg/L	0.10 mmol/L
Rheumatoid Factor ²¹	500 IU/mL	500 IU/mL
Salicylic Acid	60 mg/dL	4.34 mmol/L
Sulfasalazine	40 µg/mL	100 µmol/L
Theophylline	4 mg/dL	222 µmol/L
Triamterene	60 µg/mL	237 µmol/L
Triglyceride ²²	1266 mg/dL	14 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

m. Interference testing for Dextran 40 was performed only at a TBL concentration of 0.9 mg/dL [15 µmol/L].

n. Testing performed at a TBL concentration of 25.0 ± 2.0 mg/dL [428 ± 34 µmol/L].

o. Testing performed only at a TBL concentration of 1.4 mg/dL [24 µmol/L].

Analytical Sensitivity: 0.1 mg/dL [2 µmol/L]


The analytical sensitivity represents the lower limit of detection (lowest concentration of total bilirubin that can be distinguished from zero). This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the low level (0 mg/dL [0 µmol/L]) TBL/DBI Calibrator. Purified Water Diluent (Cat. No. 710615901) or Reagent grade water should be used as the Level 1 calibrator.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics Inc.

©2008 Siemens Healthcare Diagnostics Inc.
All rights reserved.

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