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REF DF27

PI-LAB-005/00(01/10/2560)

**Dimension<sup>®</sup>** clinical chemistry system

**Flex<sup>®</sup>** reagent cartridge

**CHOL**

See shaded sections| Updated information from 2004-06 version.

Issue Date 2008-03-04


### Cholesterol

**Intended Use:** The CHOL method used on the Dimension<sup>®</sup> clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of total cholesterol in human serum and plasma.

**Summary:** The CHOL method is based on the principle first described by Stadtman<sup>1</sup> and later adapted by other workers,<sup>2,3</sup> including Rautela and Liedtke.<sup>4</sup>

Lipids and lipoproteins in circulation have been strongly associated with coronary heart disease (CHD), associated lipid metabolism disorders, and atherosclerosis, a cause of CHD.<sup>5</sup>

**Principles of Procedure:** Cholesterol esterase (CE) catalyzes the hydrolysis of cholesterol esters to produce free cholesterol which, along with preexisting free cholesterol, is oxidized in a reaction catalyzed by cholesterol oxidase (CO) to form cholest-4-ene-3-one and hydrogen peroxide. In the presence of horseradish peroxidase (HPO), the hydrogen peroxide thus formed is used to oxidize N,N-diethylaniline-HCl/4-aminoantipyrine (DEA-HCl/AAP) to produce a chromophore that absorbs at 540 nm. The absorbance due to oxidized DEA-HCl/AAP is directly proportional to the total cholesterol concentration and is measured using a polychromatic (452, 540, 700 nm) endpoint technique.

### Reagents

Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>	Source
1 – 3	Tablet <sup>c</sup>	CE	0.7 U/mL	Fungal
		CO	0.1 U/mL	Microbial
		HPO	2.4 U/mL	Plant
1 – 3	Tablet <sup>c</sup>	AAP	4.5 µmol	
		Buffer		
4 – 6	Liquid	Cholate		
		DEA	5.8 µmol	
		Surfactant		

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

c. Tablet contains excipients.

**Precautions:** Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

**Reagent Preparation:** Hydrating, diluting and mixing are automatically performed by the instrument.

**Store at:** 2 – 8 °C

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

**Open Well Stability:** 5 days for wells 1 – 6

**Specimen Collection and Handling:** Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method.<sup>5</sup>

Corvac<sup>®</sup> and SST<sup>®</sup> collection tubes do not affect the CHOL method.

Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.<sup>6</sup>

Complete clot formation should take place before centrifugation.<sup>7</sup>

Follow the instructions provided with your specimen collection device for use and processing.<sup>8</sup>

Specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.<sup>7</sup> Repeated freezing and thawing should be avoided.

Corvac<sup>®</sup> is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.

SST<sup>®</sup> is a registered trademark of Becton-Dickinson, Rutherford, NJ.

### Procedure

Materials Provided
CHOL Flex <sup>®</sup> reagent cartridge, Cat. No. DF27
Materials Required But Not Provided
CHOL Calibrator, Cat. No. DC16
Quality Control Materials

### Test Steps

Sampling,<sup>9</sup> reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension<sup>®</sup> system. For details of this processing, refer to your Dimension<sup>®</sup> Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

### Test Conditions

Sample Size	3 µL
Reagent 1 Volume	88 µL
Reagent 2 Volume	26 µL
Diluent Volume	241 µL
Temperature	37 °C
Wavelength	452, 540 and 700 nm
Type of Measurement	Polychromatic endpoint

### Calibration

Assay Range	50 – 600 mg/dL [1.3 – 15.5 mmol/L] <sup>e</sup>
Calibration Material	CHOL Calibrator (Cat. No. DC16)
Calibration Scheme	3 levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.0259) = [mmol/L]
Typical Calibration Levels	50, 250, 450 mg/dL [1.3, 6.5, 11.6 mmol/L]
Calibration Frequency	Every 3 months for any one lot Every new reagent cartridge lot
A new calibration is required	<ul style="list-style-type: none"> <li>For each new lot of Flex<sup>®</sup> reagent cartridges</li> <li>After major maintenance or service, if indicated by quality control results</li> <li>As indicated in laboratory quality control procedures</li> <li>When required by government regulations</li> </ul>
Assigned Coefficients	C <sub>1</sub> 0.3162 C <sub>2</sub> 0.7886

e. Système International d'Unités [SI Units] are in brackets.

### Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known cholesterol concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

**Results:** The instrument automatically calculates and prints the concentration of cholesterol in mg/dL [mmol/L] using the calculation scheme illustrated in your Dimension<sup>®</sup> Operator's Guide.

**Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.**

**Analytical Measurement Range (AMR):** 50 – 600 mg/dL [1.3 – 15.5 mmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

• Samples with results in excess of 600 mg/dL [15.5 mmol/L] should be repeated on dilution.

**Manual Dilution:** Make appropriate dilutions with Reagent grade water to obtain result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

**Autodilution (AD):** Refer to your Dimension<sup>®</sup> Operator's Guide.

CHOL results less than 50 mg/dL [1.3 mmol/L] should be reported as "less than 50 mg/dL [1.3 mmol/L]" instead of the numerical value.

### Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension<sup>®</sup> Operator's Guide.

For linearity verification use high and low serum pools or commercially available linearity verification products. It is recommended that the specimens used to determine the linearity have triglycerides less than 600 mg/dL [6.78 mmol/L].

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
220 mg/dL [5.7 mmol/L]	>5 mg/dL [0.13 mmol/L]
400 mg/dL [10.3 mmol/L]	>7 mg/dL [0.18 mmol/L]

### Interfering Substances

Potassium Oxalate/Sodium Fluoride can decrease cholesterol results an average of 12%.

Li Heparin can depress cholesterol results by an average of 4 mg/dL [0.1 mmol/L] at a level of 200 mg/dL [5.2 mmol/L].

Bilirubin (conjugated) of 8.1 mg/dL [139 µmol/L] and bilirubin (unconjugated) of 9.4 mg/dL [161 µmol/L] decrease the CHOL result by 15 mg/dL [0.4 mmol/L] at CHOL concentration of 150 mg/dL [3.9 mmol/L].

Bilirubin (conjugated) of 12.8 mg/dL [219 µmol/L] and bilirubin (unconjugated) of 14.7 mg/dL [251 µmol/L] decrease the CHOL result by 25 mg/dL [0.7 mmol/L] at CHOL concentration of 250 mg/dL [6.5 mmol/L].

Bilirubin (unconjugated) of 20 mg/dL [342 µmol/L] decreases a CHOL result of 178 mg/dL [4.6 mmol/L] by 15%.

Hemoglobin (hemolysate) of 1000 mg/dL [0.62 mmol/L] (monomer) decreases a CHOL result of 177 mg/dL [4.6 mmol/L] by 15%.

Lipemia (Intralipid<sup>®</sup>) at 1000 mg/dL [11.3 mmol/L] tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid<sup>®</sup> is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

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#### Specific Performance Characteristics<sup>f</sup>

##### Precision<sup>a</sup>

Material	Mean mg/dL [mmol/L]	Standard Deviation (% CV)	
		Within-run	Between-day
QCS™ Control			
Serum Assayed			
Normal	143 [3.7]	1.2 [0.03] (0.84)	1.9 [0.05] (1.3)
Abnormal	105 [2.7]	0.97 [0.03] (0.92)	1.8 [0.05] (1.8)

- f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).  
g. Specimens at each level were analyzed in triplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

QCS™ is a trademark of Gifford, Irvine, CA 92714.

##### Method Comparison

##### Regression Statistics<sup>a</sup>

Comparative Method	Slope	Intercept mg/dL [mmol/L]	Correlation Coefficient	n
CHOL method on the aca® discrete clinical analyzer	1.038	-3.6 [-0.1]	0.995	86

h. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept

#### Specificity

##### HIL Interference

The CHOL method was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	CHOL Concentration mg/dL [mmol/L]	Bias <sup>i</sup> %
Hemoglobin (hemolysate)	500 mg/dL [0.31 mmol/L] (monomer)	177 [4.6]	<10
Bilirubin (unconjugated)	5 mg/dL [86 µmol/L]	178 [4.6]	<10
Lipemia (Intralipid®)	600 mg/dL [6.86 mmol/L]	175 [4.5]	<10

i. Analyte results should not be corrected based on this bias

#### Non-Interfering Substances

The following substances have a negligible effect (<10%) on a CHOL result of 217 mg/dL [222.6 mmol/L] at the concentrations indicated.

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	17.6 µmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	5 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	12 mg/dL	251 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lipemia	600 mg/dL	6.78 mmol/L (endogenous triglyceride)
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	421 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Interference from the statins at the concentrations indicated, when added to a 151 mg/dL [3.9 mmol/L] CHOL serum pool respectively, is less than 10%.

Substance	Test Concentration	SI Units
Atorvastatin	3.6 µg/mL	2.98 µmol/L
Lovastatin	2.0 µg/mL	4.94 µmol/L
Pravastatin	10.3 µg/mL	23.07 µmol/L
Simvastatin	6.9 µg/mL	16.48 µmol/L

##### Certification

This method has been evaluated by and met the certification acceptance criteria of the Cholesterol Reference Method Laboratory Network (CRMLN).

**Analytical Sensitivity: 50 mg/dL [1.3 mmol/L]**

The analytical sensitivity represents the lowest concentration of CHOL that can be distinguished from zero. This sensitivity is defined as the low assay range.

**Symbols Key:** See adjacent panel.

**Bibliography:** See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics Inc.

©2008 Siemens Healthcare Diagnostics Inc.  
All rights reserved.

[illegible]