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REF DF48B

PI-LAB-007/00(01/10/2560)

Dimension® clinical chemistry system

Flex® reagent cartridge

AHDL

Issue Date 2009-07-22

Automated HDL Cholesterol

Intended Use: The AHDL method is an *in vitro* diagnostic test for the quantitative measurement of high-density lipoprotein cholesterol (HDL-C) in human serum and plasma on the Dimension® clinical chemistry system. Measurements of HDL-C are used as an aid in the diagnosis of lipid disorders (such as diabetes mellitus), various liver and renal diseases and in the assessment of risk for atherosclerosis and cardiovascular disease.


Summary: Plasma lipoproteins are spherical particles of varying composition. The outer surface of these particles is made up of phospholipids, free cholesterol and protein; the inner core contains mostly esterified cholesterol and triglyceride. Lipoproteins function to solubilize and transport cholesterol and triglycerides in the blood stream.

Four types of lipoproteins are recognized clinically based on the relative proportions of their lipid and protein content: chylomicrons, very low-density lipoproteins (VLDL), low-density lipoproteins (LDL) and high-density lipoproteins (HDL). The primary function of HDL is to transport cholesterol from peripheral tissues to the liver where it is metabolized.¹ This process, known as reverse cholesterol transport, has been proposed to be a cardiovascular protective mechanism. Patients with low levels of HDL cholesterol are generally considered to be at increased risk for coronary artery disease.²⁻⁴

The determination of serum HDL cholesterol level is a useful tool in identifying at-risk patients. The Adult Treatment Panel of the National Cholesterol Education Program (NCEP) recommends that all adults 20 years of age and over should have their total and HDL cholesterol levels measured at least every 5 years to screen for risk of coronary artery disease.⁵

The reference method for measuring HDL cholesterol utilizes ultracentrifugation and chemical precipitation to separate HDL from the other lipoproteins, followed by a determination of the cholesterol content using the Abell-Kendall method. Since this procedure is tedious and requires an ultracentrifuge, most laboratories now use methods which selectively measure HDL-C by chemical means.

Principles of Procedure: The AHDL assay measures HDL cholesterol levels directly without the need for sample pretreatment or specialized centrifugation steps, using a two reagent format. In the first reaction, chylomicrons, VLDL and LDL form water soluble complexes with dextran sulfate in the presence of magnesium sulfate. These complexes are resistant to the polyethylene glycol (PEG)-modified cholesterol esterase and cholesterol oxidase that react with HDL cholesterol. In the presence of oxygen, the HDL cholesterol is oxidized to Δ^4 -cholestenone and hydrogen peroxide. The generated hydrogen peroxide then reacts with 4-aminoantipyrine and sodium N-(2-hydroxy-3-sulfoethyl)-3,5-dimethoxyaniline (HSDA) in the presence of peroxidase to form a colored dye that is measured using a bichromatic (600/700 nm) endpoint technique. The color intensity of the dye is directly proportional to the serum HDL-C concentration.

Reagents				
Wells ^a	Form	Ingredient	Concentration ^b	Source
1, 2, 3 (Reagent 1)	Liquid	HEPES Buffer	10.07 mmol/L, pH 7.4	
		2-(N-cyclohexylamino)-ethanesulfonic acid	96.95 mmol/L	
		Dextran Sulfate	1.5 g/L	
		Magnesium Nitrate Hexahydrate	≥ 11.7 mmol/L	
		N-(2-hydroxy-3-sulfoethyl)-3,5-dimethoxyaniline	0.96 mmol/L	
		Ascorbate Oxidase	≥ 50 µkat/L	bacterial
		Peroxidase	≥ 16.7 µkat/L	horseradish
		Preservative		
		HEPES Buffer	10.07 mmol/L, pH 7.0	
		PEG-Cholesterol Esterase	≥ 3.33 µkat/L	bacterial
4 (Reagent 2)	Liquid	PEG-Cholesterol Oxidase	≥ 127 µkat/L	bacterial
		Peroxidase	≥ 333 µkat/L	horseradish
		4-Aminoantipyrine	2.46 mmol/L	
		Preservative		
5, 6	Liquid	NaOH ^c	1.00 M	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per well in a cartridge.

c. Sodium hydroxide is used as a probe cleaning solution and is not used in the reaction.

Risk and Safety

Corrosive. Contains Sodium Hydroxide.
R34: Causes burns.
S26: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

S36/37/39: Wear suitable protective clothing, gloves and eye/face protection.
S45: In case of accident or if you feel unwell, seek medical advice immediately. (Show the label where possible.)

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Contains sodium azide (< 0.1%) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations. Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 3
10 days for well 4
15 days for wells 5 – 6

Specimen Collection and Handling

Recommended specimen types: serum and plasma (lithium or sodium heparin).

Serum and plasma should be collected after a 12-hour period of fasting using recommended procedures for collection of diagnostic blood specimens by venipuncture.^{6,7} Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.⁸ Specimens should be free of particulate matter.

Serum or plasma samples may be refrigerated at 2 – 8 °C for up to 7 days if not tested within 8 hours. For longer storage, samples may be frozen at -70 °C for up to three months. Samples may be frozen once. When using samples stored at 2 – 8 °C or -70 °C, increases or decreases of up to 10% may be observed in HDL-C concentrations.

Follow the instructions provided with your specimen collection device for use and processing.⁹

Procedure

Materials Provided

AHDL Flex® reagent cartridge, Cat. No. DF48B

Materials Required But Not Provided

AHDL Calibrator, Cat. No. DC48B
Quality Control Materials

Test Steps

Sampling⁹, reagent delivery, mixing, and processing are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

Test Conditions

Sample Volume	3 µL
Reagent 1 Volume	300 µL
Reagent 2 Volume	100 µL
Temperature	37.0 °C ± 0.1 °C
Reaction time	8.6 minutes
Wavelength	600 and 700 nm
Type of Measurement	Bichromatic end point

Calibration

The general calibration procedure is described in your Dimension® Operator's Guide. The following information should be considered when calibrating the AHDL method:

Measurement Range	3 – 150 mg/dL [0.08 – 3.89 mmol/L] ^e
Calibration Material	AHDL Calibrator, Cat. No. DC48B
Calibration Scheme	3 levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.0259) = [mmol/L]

Typical Calibration Levels
0, 40, 165 mg/dL
[0.00, 1.04, 4.27 mmol/L]

Calibration Frequency
Every 90 days for any one lot

A new calibration is required

- For each new lot of Flex® reagent cartridges
- After major maintenance or service, if indicated by quality control results
- As indicated in laboratory quality control procedures
- When required by government regulations

Assigned Coefficients	C ₀ 0.000
	C ₁ 1.700

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Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known HDL-C concentrations. The National Cholesterol Education Program (NCEP) Lipid Standardization Panel (LSP) recommends two levels of controls: one in the normal range (40 – 60 mg/dL [1.04 – 1.55 mmol/L]) and one near the concentrations for decision-making (<40 mg/dL [<1.04 mmol/L]).¹ Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of AHDL in mg/dL [mmol/L] using the calculation scheme described in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range: 3 – 150 mg/dL [0.08 – 3.89 mmol/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

- Samples with results in excess of 150 mg/dL [3.89 mmol/L] are reported as "Above Assay Range" and should be repeated on dilution.

Manual Dilution: Dilute with Reagent Grade water to obtain results within reportable range. Recommended dilution is 1:2. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): The recommended autodilute sample volume is 2 µL. Refer to your Dimension® Operator's Guide.

- Samples with results less than 3 mg/dL [0.08 mmol/L] should be reported as "less than 3 mg/dL [0.08 mmol/L]".

Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in AHDL results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual and not reported.

A system malfunction may exist if the following five-test precision is observed:

AHDL Concentration		SD
45 mg/dL [1.17 mmol/L]		>1.4 mg/dL [0.04 mmol/L]
60 mg/dL [1.55 mmol/L]		>2.0 mg/dL [0.05 mmol/L]

Interfering Substances

Abnormal liver function may affect lipid metabolism; therefore, HDL cholesterol and LDL cholesterol results are of limited diagnostic value. In some patients with abnormal liver function, the HDL-C result may significantly differ from the designated comparison method (DCM) result.¹

In very rare cases, gammopathy, in particular type IgM (Waldenstrom's macroglobulinemia), may cause unreliable results.^{10,11}

Expected Values

The National Cholesterol Education Program Adult Treatment Panel III (NCEP-ATP III) provides the following classifications of HDL-C concentrations:

HDL < 40 mg/dL [1.04 mmol/L] Low HDL Cholesterol
HDL ≥ 60 mg/dL [1.55 mmol/L] High HDL Cholesterol

Each laboratory should establish its own expected values for AHDL as performed on the Dimension® system.

Specific Performance Characteristics

The following data represent typical performance for the Dimension® System.

Material	Precision ^{12,f}		
	Mean mg/dL [mmol/L]	Standard Deviation (% CV)	
		Repeatability	Within-Lab
Bio-Rad Liquid Unassayed Multiquel®			
Level 1	25 [0.65]	0.6 [0.02] (2.3)	0.6 [0.02] (2.4)
Level 2	44 [1.14]	0.6 [0.02] (1.4)	0.7 [0.02] (1.7)
Level 3	58 [1.50]	1.2 [0.03] (2.2)	1.3 [0.03] (2.2)
Serum Pool	48 [1.24]	0.9 [0.02] (1.9)	1.0 [0.03] (2.1)
Serum Pool	92 [2.38]	1.3 [0.03] (1.4)	1.7 [0.04] (1.8)

f. CLSI EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Multiquel® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618.

Method Comparison¹³ Regression Statistics^h

Comparative Method	Slope	Intercept mg/dL [mmol/L]	Correlation Coefficient	n
Dimension® AHDL method DF48A	1.04	-3.9 [-0.10]	1.00	130 ^g

g. CLSI EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.

h. The range of AHDL values in the correlation study was 4 – 131 mg/dL [0.10 – 3.39 mmol/L].

Specificity

HIL Interference

The AHDL method was evaluated for interference according to CLSI EP7-A2.¹⁴

Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	HDL Cholesterol mg/dL [mmol/L]	Bias* %
Hemoglobin (hemolysate)	Hemoglobin (monomer) 1000 mg/dL [0.62 mmol/L]	37 [0.96]	< 10
Bilirubin (unconjugated)	80 mg/dL [1368 µmol/L]	43 [1.11]	< 10
Bilirubin (conjugated)	60 mg/dL [1026 µmol/L]	43 [1.11]	< 10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	37 [0.96]	< 10

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

* **Analyte results should not be corrected based on this bias.**

Non Interfering Substances

The following substances do not interfere with the AHDL method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at HDL-C concentrations of 40 mg/dL [1.04 mmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Albumin	6000 mg/dL	60 g/L
Amikacin	8.0 mg/dL	137 µmol/L
Ampicillin	5 mg/dL	143 µmol/L
Ascorbic Acid	6.0 mg/dL	342 µmol/L
Atorvastatin	3.6 µg/mL	2.98 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Clofibrate	400 µg/mL	1.65 mmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.51 mg/dL	17.9 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Fenofibrate	46 µg/mL	127.4 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gemfibrozil	124 µg/mL	495.3 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	6 g/dL	60 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Lovastatin	2 µg/mL	4.94 µmol/L
Niacin	31 µg/mL	251.8 µmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Pravastatin	10.3 µg/mL	23.07 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Rheumatoid Factor	775 IU/mL	775 IU/mL
Salicylic Acid	60 mg/dL	4.34 mmol/L
Simvastatin	6.9 µg/mL	16.48 µmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	1150 mg/dL	13.0 mmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

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Cross-reactivity

Low-density lipoprotein cholesterol was evaluated for cross-reactivity in the presence of 40 mg/dL [1.04 mmol/L] HDL-C. The percent cross-reactivity was calculated as follows:

$$\% \text{ crossreactivity} = \frac{\text{measured HDL-C} - \text{control HDL-C}}{\text{LDL-C}} \times 100\%$$

LDL-C Concentration	% Crossreactivity
282 mg/dL	0.3

Limit of Detection¹: 3 mg/dL [0.08 mmol/L]

The limit of detection represents the lowest concentration of high density lipoprotein cholesterol that can be detected with at least 95% probability. CLSI/NCCLS EP17-A, Protocols for Determination of Limits of Detection and Limits of Quantitation, was followed using 60 replicate determinations at two concentration levels to determine Limit of Blank and Limit of Detection.

Certification

This method has been evaluated by and met the certification acceptance criteria of the Cholesterol Reference Method Laboratory Network (CRM/LN).¹

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension® and Flex® are trademarks of Siemens Healthcare Diagnostics.

©2009 Siemens Healthcare Diagnostics Inc.

All rights reserved.

FM-LAB-008/00(01/10/2560)