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#### Specific Performance Characteristics<sup>a</sup>

Material	Precision <sup>d</sup>		Standard Deviation (% CV)	
	Mean mg/dL [mmol/L]	Within-run	Total	
QCS™ Control Serum Assayed				
Normal	2.2 [0.91]	0.07 [0.03] (3.2)	0.11 [0.04] (5.0)	
Abnormal	3.9 [1.60]	0.18 [0.07] (4.6)	0.19 [0.08] (4.9)	
Fisher Urine Chemistry Control				
Urine				
Level 1	20.1 [8.27]	0.4 [0.16] (2.0)	0.6 [0.25] (3.0)	
Level 2	4.2 [1.73]	0.1 [0.04] (2.4)	0.1 [0.04] (2.4)	
Urine Pool	9.8 [4.03]	0.1 [0.04] (1.0)	0.2 [0.08] (2.0)	

e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (see Dimension® Operator's Guide).

f. Specimens at each level were analyzed in duplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

QCS™ is a trademark of Gilford, Irvine, CA 92714.

#### Method Comparison

Comparative Method	Regression Statistics <sup>a</sup>		Correlation Coefficient	n
	Slope	Intercept mg/dL [mmol/L]		
Serum				
aca® discrete clinical analyzer	0.98	-0.14 [-0.06]	0.970	110
Urine				
Johnson and Johnson Vitros clinical analyzer	1.00	-0.11 [-0.04]	0.987	58
aca® discrete clinical analyzer	1.016	-0.04 [-0.02]	0.990	109

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#### Specificity

##### HIL Interference

The MG method was evaluated for interference from hemolysis, icterus, and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration SI Units	MG Concentration mg/dL [mmol/L]	Bias <sup>h</sup> %
Hemoglobin (hemolysate)	200 mg/dL [0.12 mmol/L] (monomer)	0.94 [0.39]	<10
Bilirubin (unconjugated)	40 mg/dL [684 µmol/L]	0.97 [0.40]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	0.92 [0.38]	<10

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REF DF57

PI-LAB-013/00(01/10/2560)

**Dimension®** clinical chemistry system

**Flex®** reagent cartridge

**MG**

See shaded sections: Updated information from 2003-07 version.


Issue Date 2008-10-02

### Magnesium

**Intended Use:** The MG method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of magnesium in human serum, heparinized plasma and urine.

**Summary:** The magnesium method is a modification of the methylthymol blue (MTB) complexometric procedure described by Connerty, Lau, and Briggs.<sup>1</sup> The barium salt of ethylenedis (oxyethylenetri) tetraacetic acid (Ba-EGTA) is used to reduce interference due to calcium which also reacts with MTB.<sup>2</sup>

**Principles of Procedure:** MTB forms a blue complex with magnesium. Calcium interference is minimized by forming a complex between calcium and Ba-EGTA (chelating agent). The amount of MG-MTB complex formed is proportional to the magnesium concentration and is measured using a bichromatic (600 and 510 nm) endpoint technique. For Dimension® AR (software version 4.4 and above), ES (software 2.7 or 4.3 and above) and XL/RxL (software version 5.0 and above), a sample blank is used to minimize bilirubin interference.

#### Reagents

Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>
1 – 3	Liquid	MTB Acetic acid Potassium Sorbate	0.0528 g/L
4 – 6	Liquid	Ba-EGTA Sodium metaborate Buffer Microbial inhibitors	0.5 mM

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

**Precautions:** Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use

**Reagent Preparation:** All reagents are liquid and ready to use.

**Store at:** 2 – 8 °C

**Expiration:** Refer to shelf carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

**Open Well Stability:** 2 days for wells 1 – 8

**Specimen Collection and Handling:** Normal procedures for collecting serum, plasma and urine may be used for samples to be analyzed by this method.<sup>3,4</sup>

Follow the instructions provided with your specimen collection device for use and processing.<sup>5</sup>

Complete clot formation should take place before centrifugation.<sup>6</sup>

Separated specimens are stable for 3 days at room temperature, 3 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder for up to a year.<sup>7</sup>

A 24-hour urine specimen should be collected in a bottle containing 10 mL of 12 M hydrochloric acid.<sup>7</sup> If preservatives are not used during collection, acidify the sample to below pH 3.0 before analysis.<sup>8</sup> Acidified urine samples are stable at 4 – 6 °C for 1 week.<sup>7</sup>

Corvac® and SST® collection tubes do not affect the MG method.

Blood collection tubes containing lithium heparin, potassium oxalate or sodium fluoride do not interfere with this MG method.

EDTA of 200 mg/dL [2 g/L]<sup>c</sup> decreases the MG result by 0.4 mg/dL [0.16 mmol/L] at a magnesium concentration of 1.8 mg/dL [0.74 mmol/L].

c. Système International d'Unités [SI Units] are in brackets.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO. SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

#### Procedure

Materials Provided
MG Flex® reagent cartridge, Cat. No. DF57
Materials Required But Not Provided
CHEM II Calibrator, Cat. No. DC20 Quality Control Materials

#### Test Steps

Sampling,<sup>4</sup> reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

#### Test Conditions

Sample Size	4 µL
Reagent 1 Volume	100 µL
Reagent 2 Volume	200 µL
Diluent Volume	696 µL
Temperature	37 °C
Wavelength	600 and 510 nm
Type of Measurement	Bichromatic endpoint

#### Calibration

Assay Range	0.0 – 20.0 mg/dL [0.0 – 8.22 mmol/L]
Calibration Material	CHEM II Calibrator, Cat. No. DC20
Calibration Scheme	3 levels, n = 3
Units	mg/dL [mmol/L] (mg/dL x 0.4114) = [mmol/L]
Typical Calibration Levels	1.0, 9.0, 18.0 mg/dL [0.41, 3.70, 7.41 mmol/L]
Calibration Frequency	Every 3 months for any one lot
A new calibration is required	<ul style="list-style-type: none"> <li>For each new lot of Flex® reagent cartridges</li> <li>After major maintenance or service, if indicated by quality control results</li> <li>As indicated in laboratory quality control procedures</li> <li>When required by government regulations</li> </ul>
Assigned Coefficients	C <sub>0</sub> -0.200 C <sub>1</sub> 0.100

#### Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known magnesium concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

**Results:** The instrument automatically calculates and prints the concentration of magnesium in mg/dL [mmol/L] using the calculation scheme illustrated in the Principles of Operation section of the Dimension® Operator's Guide.

**Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.**

**Analytical Measurement Range (AMR):** 0.0 – 20.0 mg/dL [0.00 – 8.22 mmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 20.0 mg/dL [8.22 mmol/L] should be repeated on dilution.

<b>Manual Dilution:</b>	Dilute with Reagent grade water or equivalent to produce a sample concentration within the reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.
<b>Urine:</b>	Dilute with Reagent grade water to obtain results within the reportable range. Enter dilution factor on the instrument. Reassay. Resulting readout is corrected for dilution.
<b>Autodilute (AD):</b>	Not available for this method.

#### Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to the Dimension® Operator's Guide for troubleshooting.

A system malfunction may exist if the following 5-test precision is observed:

Concentration	SD
1.0 mg/dL [0.41 mmol/L]	>0.11 mg/dL [0.05 mmol/L]
9.0 mg/dL [3.7 mmol/L]	>0.56 mg/dL [0.23 mmol/L]

#### Interfering Substances

Because MG concentration is 3 times higher in red blood cells than in serum, hemolyzed samples may give spuriously elevated magnesium results. Bias from hemolysis may vary due to individual sample variations in intracellular magnesium.<sup>9</sup>

Hemoglobin (hemolysate) of 300 mg/dL [0.19 mmol/L] (monomer) increases a magnesium concentration of 0.94 mg/dL [0.39 mmol/L] by 18%.

Bilirubin (unconjugated) of 60 mg/dL [1026 µmol/L] increases a magnesium concentration of 0.97 mg/dL [0.40 mmol/L] by 14%.

Lipemia (Intralipid®) of 3000 mg/dL [33.9 mmol/L] increases a magnesium concentration of 0.92 mg/dL [0.38 mmol/L] by 14%.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

<b>Expected Values:</b>	<b>Serum: 1.8 – 2.4 mg/dL [0.74 – 0.99 mmol/L]</b> <b>Urine: 24 – 255 mg/24hr [0.99 – 10.5 mmol/24hr]<sup>8</sup></b>
-------------------------	--

The reference interval was calculated non-parametrically and represents the central 95% of results determined from a population of healthy adults (n = 118; 61 males, ages 21 – 56 and 57 females, ages 20 – 55).

Each laboratory should establish its own reference interval for magnesium as performed on the Dimension® system.

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#### Specific Performance Characteristics\*

Material	Precision <sup>†</sup>		Standard Deviation (% CV)	
	Mean mg/dL [mmol/L]	Within-run	Total	
QCS™ Control Serum Assayed				
Serum				
Normal	2.2 [0.91]	0.07 [0.03] (3.2)	0.11 [0.04] (5.0)	
Abnormal	3.9 [1.60]	0.18 [0.07] (4.6)	0.19 [0.08] (4.9)	
Fisher Urine Chemistry Control				
Urine				
Level 1	20.1 [8.27]	0.4 [0.16] (2.0)	0.6 [0.25] (3.0)	
Level 2	4.2 [1.73]	0.1 [0.04] (2.4)	0.1 [0.04] (2.4)	
Urine Pool	9.8 [4.03]	0.1 [0.04] (1.0)	0.2 [0.08] (2.0)	

- e. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (see Dimension® Operator's Guide).
- f. Specimens at each level were analyzed in duplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

QCS™ is a trademark of Gifford, Irvine, CA 92714.

Method Comparison				
Comparative Method	Regression Statistics <sup>d</sup>		Correlation Coefficient	n
	Slope	Intercept mg/dL [mmol/L]		
Serum				
aca® discrete clinical analyzer	0.98	-0.14 [-0.06]	0.970	110
Serum				
Johnson and Johnson Vitros clinical analyzer	1.00	-0.11 [-0.04]	0.987	58
Urine				
aca® discrete clinical analyzer	1.016	-0.04 [-0.02]	0.990	109

g. Model equation for regression statistics is: [results of Dimension® system] = slope x [comparative method results] + intercept.

#### Specificity

HIL Interference			
The MG method was evaluated for interference from hemolysis, icterus, and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".			
Substance Tested	Test Concentration SI Units	MG Concentration mg/dL [mmol/L]	Bias <sup>†</sup> %
Hemoglobin (hemolysate)	200 mg/dL [0.12 mmol/L] (monomer)	0.94 [0.39]	<10
Bilirubin (unconjugated)	40 mg/dL [684 µmol/L]	0.97 [0.40]	<10
Lipemia (Intralipid®)	1000 mg/dL [11.3 mmol/L]	0.92 [0.38]	<10

h. Analyte results should not be corrected based on this bias.

#### Non-Interfering Substances

The following substances do not interfere with the MG method when present in plasma in the amounts indicated. Systematic inaccuracies (bias) due to these substances are less than 10% at following concentrations of magnesium tested.

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Calcium	20 mg/dL	4.99 mmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	252 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	20 µg/dL	70 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	800 mg/dL	174 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	16 µg/mL	29.4 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	25U/mL	25000 U/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Pencillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	80 mg/mL	344 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	100 µg/dL	555 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L
Uric Acid	20 mg/dL	1.2 mmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

#### Analytical Sensitivity: 0.0 mg/dL [0.00 mmol/L]


The analytical sensitivity represents the lowest concentration of magnesium that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water 0.0 mg/dL [0.00 mmol/L].

**Symbols Key:** See adjacent panel.

**Bibliography:** See adjacent panel.

Dimension®, aca® and Flex® are trademarks of Siemens Healthcare Diagnostics Inc.

©2008 Siemens Healthcare Diagnostics Inc.  
All rights reserved.

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