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a. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed on Dimension® (Refer to your Dimension® Operator's Guide).

f. Reproducibility testing was done in accordance with the CLSI/NCCLS Approved Guideline for Evaluation of Precision Performance of Clinical Chemistry Devices (EP5-A, Feb. 1999).

g. Specimens at each level were analyzed in duplicate, once a day, for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

HDL Plus® is a registered trademark of Quantimetrix Corporation, Redondo Beach, CA.

Method Comparison

Regression Statistics^a

Comparative Method	Slope	Intercept mg/dL (mmol/L)	Correlation Coefficient	n ⁱ
N-geneous™ LDL Cholesterol Reagent ^j	0.95	4.7 [0.12]	0.997	122
Beta-Quantification ^k	1.01	3.3 [0.08]	0.982	49

h. Model equation for regression statistics is: Result of Dimension® analyzer = (Slope x Comparative method result) + Intercept.

i. The range of LDL-C values in the correlation study was 35 to 273 mg/dL [0.90 to 7.06 mmol/L] for the N-geneous™ Cholesterol Reagent^j and 70 to 246 mg/dL [1.81 to 6.36 mmol/L] for Beta-Quantification. Correlation study was performed using serum samples.

j. N-geneous™ LDL Cholesterol Reagent^j is manufactured by Genzyme Corporation, Cambridge, MA 02139-1562. Testing was performed on the Beckman CX9 chemistry analyzer manufactured by Beckman Coulter Inc., Brea, CA 92621.

k. Beta-Quantification Reference Method for LDL-C is a three-step procedure^F involving ultracentrifugation, precipitation with heparin-manganese reagent and quantification with the Abel-Kendall reference method for cholesterol. This study was performed by the Core Laboratory for Clinical Studies at Washington University, St. Louis, Missouri, 63110.

This method has been certified using the Low Density Lipoprotein (LDL) Cholesterol Method Evaluation Protocol by the Cholesterol Reference Method Laboratory Network (CRMLN) at Northwest Lipid Research Laboratories, University of Washington, Seattle, Washington, 98103.

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REF DF41A

PI-LAB-018/00 (01/10/2560)

Dimension® clinical chemistry system

Flex® reagent cartridge

AST

See shaded sections; Updated information from 2003-07 version.


Issue Date 2008-02-12

Aspartate Aminotransferase

Intended Use: The AST method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of aspartate aminotransferase activity in human serum or plasma.

Summary: The aspartate aminotransferase method is an adaptation of the methodology recommended by the International Federation of Clinical Chemistry (IFCC).¹ The method uses the coenzyme pyridoxal-5-phosphate (P5P) to activate the apoenzyme and lactic acid dehydrogenase (LDH) to eliminate pyruvate interference. Significant elevations of AST are found in diseases of the liver such as hepatitis, necrosis, jaundice, and cirrhosis. AST levels can be elevated even before clinical jaundice appears.²

Principles of Procedure: Aspartate aminotransferase (AST) catalyzes the transamination from L-aspartate to α-ketoglutarate, forming L-glutamate and oxalacetate. The oxalacetate formed is reduced to malate by malate dehydrogenase (MDH) with simultaneous oxidation of reduced nicotinamide adenine dinucleotide (NADH). The change in absorbance with time due to the conversion of NADH to NAD is directly proportional to the AST activity and is measured using a bichromatic (340, 700 nm) rate technique.

Reagents

Wells ^a	Form ^b	Ingredient	Concentration ^c	Source
1 – 3	Tablet	MDH LDH NADH	3000 U/L 900 U/L 0.21 mmol/L	porcine muscle
1 – 3	Tablet	P5P Buffer	0.18 mmol/L	
4 – 6	Tablet (2/well)	Aspartic acid α-ketoglutaric acid	180 mmol/L 11.5 mmol/L	

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Tablet contains excipients.

c. Nominal value per well when hydrated.

Risk and Safety:

Harmful. Contains 3, 5, 7-triaza-1-azoniacyclo[3.3.1.1^{3,7}] decane, 1-(3-chloro-2-propenyl)-Chloride, (z).
R22: Harmful if swallowed.
Xn S46: If swallowed, seek medical advice immediately and show this container or label.

Safety data sheets (MSDS/SDS) available on www.siemens.com/diagnostics

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: Hydrating, diluting and mixing are automatically performed by the instrument.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed or unhydrated cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 6

Specimen Collection and Handling: Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method. In the preparation of serum or plasma samples, avoid prolonged contact with the separated red cells.²

Follow the instructions provided with your specimen collection device for use and processing.³

Complete clot formation should take place before centrifugation.

Specimens are stable for 8 hours at room temperature, 2 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder.^{3,4,5}

Corvac® and SST® collection tubes do not affect the AST results.

Corvac® is a registered trademark of Monoject Division of Sherwood Medical, St. Louis, MO.

SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided

AST Flex® reagent cartridge, Cat. No. DF41A

Materials Required But Not Provided

Enzyme Verifier, Cat. No. DC19

Quality Control Materials

Test Steps

Sampling,^{4a} reagent delivery, mixing, processing, and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

d. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

e. An alternate sample size of 20 µL can be programmed; refer to the Dimension® Operator's Guide for the use of alternate sample size.

Test Conditions

Sample Size	40 µL (20 µL) ^a
Reagent 1 Volume	100 µL
Reagent 2 Volume	65 µL
Diluent Volume	235 µL
Temperature	37 °C
Wavelength	340 and 700 nm
Type of Measurement	Bichromatic rate

Verification

Assay Range (@ 37 °C)	0 – 1000 U/L
Verification Material	Enzyme Verifier (Cat. No. DC19)
Verification Scheme	3 levels, n = 3
Units	U/L
Typical Verification Levels	50, 400, 800 U/L
Verification Frequency	Every 90 days for any one lot
A new verification is required:	<ul style="list-style-type: none"> For each lot of Flex® reagent cartridges After major maintenance or service, if indicated by quality control results As indicated in laboratory quality control procedures When required by government regulations
Assigned Coefficients	Standard sample size = 40 µL C ₁ 2.000 C ₂ -3.537 Alternate sample size = 20 µL C ₁ -2.000 C ₂ -7.040

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known aspartate aminotransferase activity.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the activity of aspartate aminotransferase in U/L using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 1000 U/L

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 1000 U/L should be repeated on dilution.

Manual Dilution: Results in excess of 1000 U/L should be repeated after diluting the sample with Enzyme Diluent (Cat. No. 790035901) or equivalent to produce a sample result within the assay range. Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

Autodilution (AD): Refer to your Dimension® Operator's Guide.

Samples with results less than 6 U/L should be reported as "less than 6 U/L".

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed at the standard sample size (40 µL):

AST Activity	SD
40 U/L	>2.5 U/L
440 U/L	>8 U/L
830 U/L	>15 U/L

Interfering Substances

Hemolysis falsely elevates AST results. Hemolyzed specimens should not be used.

Lipemia (Intralipid®) of 600 mg/dL [6.78 mmol/L]^f tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

f. Système International d'Unités [SI Units] are in brackets.

Expected Values: 15 – 37 U/L

The reference interval was calculated non-parametrically and represents the central 95% of the population. (n = 245, adults)

Each laboratory should establish its own reference interval for aspartate aminotransferase as performed on the Dimension® system.

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Specific Performance Characteristics^a

Material	Mean U/L	Precision ^b	
		Standard Deviation Within-run	Total
Multiquel [®]			
Level 1	46	1.2 (2.7)	2.4 (5.2)
Level 2	190	1.6 (0.8)	3.9 (2.1)
Moni-Trol [®] ^j			
Level 1	25	2.8 (11.5)	3.0 (12.3)
Level 2	120	2.9 (2.4)	3.7 (3.1)

g. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension[®] Operator's Guide).

h. Specimens at each level were analyzed in duplicate for 20 days. The within-run and total standard deviations were calculated by the analysis of variance method.

i. Using reduced sample size (20 µL).

Multiquel[®] is a registered trademark of BioRad Corp., Irvine, CA 92714, USA.

Moni-Trol[®] is a registered trademark of Medical Analysis Systems Inc., Camarillo, CA 93012-8058, USA.

Method Comparison

Comparative Method	Slope	Regression Statistics ⁱ		n
		Intercept U/L	Correlation Coefficient	
GOT/AST method on the aca [®] discrete clinical analyzer	0.96	-9.3	0.995	144 ^k
Reduced vs. standard sample size ^l	1.02	-4.5	0.999	63 ^m

j. Model equation for regression statistics is: Results of Dimension[®] AST system = [slope x comparative method results] + intercept.

k. Range of samples: 8 – 408 U/L.

l. Model equation for regression statistics is: Results of Dimension[®] system using reduced sample size (20 µL) = [slope x results of Dimension[®] system using standard sample size (40 µL)] + intercept.

m. Range of samples: 4 – 993 U/L.

Specificity

HIL Interference			
The AST method (using the standard sample size of 40 µL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".			
Substance Tested	Test Concentration		AST Activity
	[SI Units]	U/L	
Hemoglobin (hemolysate) ^m	50 mg/dL	53	<10
Bilirubin (unconjugated)	[0.031 mmol/L] (monomer)		
	20 mg/dL	54	<10
	[342 µmol/L]		
Lipemia (Intralipid [®])	200 mg/dL	58	<10
	[2.26 mmol/L]		

The AST method (using an alternate sample size of 20 µL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration		AST Activity
	[SI Units]	U/L	
Hemoglobin (hemolysate) ^p	50 mg/dL	42	16
Bilirubin (unconjugated)	[0.031 mmol/L] (monomer)		
	20 mg/dL	36	<10
	[342 µmol/L]	36	10
	40 mg/dL		
Lipemia (Intralipid [®])	[1026 µmol/L]		
	200 mg/dL	34	<10
	[2.26 mmol/L]		

n. Analyte results should not be corrected based on this bias.

o. Because AST activity is 10 times higher in red blood cells than in serum, hemolyzed samples are likely to give spuriously elevated AST results. Bias from hemolysis may vary due to individual sample variations in intracellular AST.

Non-Interfering Substances

The following substances have been shown to have no measurable effect on the AST result at the concentrations indicated:

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.660 µmol/L
Amikacin	15 mg/dL	256 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic Acid	5 mg/dL	227 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 ng/dL	6.27 µmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2652 µmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	20 µg/dL	70 µmol/L
Digoxin	20 ng/mL	25.6 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	800 mg/dL	174 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	16 µg/dL	29.4 µmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Nortriptyline	1000 ng/mL	3797 nmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	421 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.2 mg/dL	4.91 µmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	100 µg/dL	555 µmol/L
Urea	500 mg/dL	83.3 mmol/L
Uric Acid	20 mg/dL	1190 µmol/L
Valproic Acid	50 mg/dL	3467 µmol/L

Analytical Sensitivity: 5 U/L

The analytical sensitivity represents the lowest activity of AST that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of Reagent grade water.

Symbols Key: See adjacent panel.

Dimension[®], aca[®] and Flex[®] are trademarks of Siemens Healthcare Diagnostics Inc.

©2008 Siemens Healthcare Diagnostics Inc.

All rights reserved.

FM-LAB-008/00(01/10/2560)