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Precision^g

Material	Mean	Standard Deviation (% CV)	
	mg/dL [μmol/L]	Within-run	Between-day
Serum			
QCS™ Control Serum Assayed			
Normal	5.1 [303]	0.07 [4] (1.4)	0.07 [4] (1.4)
Abnormal	9.0 [535]	0.11 [6] (1.2)	0.12 [7] (1.3)
Moni-Trol® ^h			
Level 1	3.6 [214]	0.16 [9.5] (4.4)	0.15 [8.9] (4.1)
Level 2	9.7 [577]	0.10 [5.9] (1.1)	0.16 [9.5] (1.6)
Urine			
Fisher Urine Chemistry Control			
Level 1	21.6 [1285]	0.6 [36] (2.8)	0.8 [48] (3.7)
Level 2	27.9 [1660]	0.4 [24] (1.4)	0.7 [42] (2.5)
Urine Pool	37.5 [2230]	0.5 [30] (1.3)	0.8 [48] (2.1)
Urichem Trak ^h			
Level 1	13.6 [809]	0.72 [43] (5.3)	0.88 [52] (6.5)
Level 2 (spiked)	35.4 [2106]	0.80 [48] (2.3)	2.95 [175] (8.4)

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h. Using reduced sample size (10 μL).

QCS™ is a trademark of Gilford, Irvine, CA 92714.

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Method Comparison

Regression Statistics

Comparative Method	Slope	Intercept mg/dL [μmol/L]	Correlation Coefficient	n
Serum				
URCA method on the aca® discrete clinical analyzer ^l	0.968	0.16 [9]	0.997	101
Reduced vs. Standard sample size ^l	1.030	-0.03 [-2]	0.999	60 ^k
Urine				
URCA method on the aca® discrete clinical analyzer ^l	0.997	-0.69 [-41]	0.962	114

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Specificity

HIL Interference

The URCA method (using the standard sample size of 17 μL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	URCA Concentration mg/dL [μmol/L]	Bias ^l %
Hemoglobin (hemolysate)	1000 mg/dL	6.4 [382]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL [1368 μmol/L]	6.5 [383]	<10
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	6.3 [377]	<10

The URCA method (using the alternate sample size of 10 μL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P. Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	URCA Concentration mg/dL [μmol/L]	Bias ^l %
Hemoglobin (hemolysate)	1000 mg/dL	6.7 [397]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL [1368 μmol/L]	6.7 [397]	<10
Lipemia (Intralipid®)	600 mg/dL [6.78 mmol/L]	6.6 [377]	<10

l. Analyte results should not be corrected based on this bias.

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REF DF77

PI-LAB-004/00(01/10/2560)

Dimension® clinical chemistry system

Flex® reagent cartridge

URCA

See shaded sections: Updated information from 2008-04 version

Issue Date 2009-04-01

Uric Acid

Intended Use: The URCA method used on the Dimension® clinical chemistry system is an *in vitro* diagnostic test intended for the quantitative determination of uric acid in human serum, plasma and urine.

Summary: The uric acid method is a modification of the uricase method first reported by Bulger and Johns,¹ later modified by Kalcikar.² Measurement of uric acid by monitoring the loss of absorbance at 293 nm following uricase treatment is generally recognized as being more specific and less subject to interference than other, indirect methods.

Principles of Procedure: Uric acid, which absorbs light at 293 nm is converted by uricase to allantoin, which is nonabsorbing at 293 nm. The change in absorbance at 293 nm due to the disappearance of uric acid is directly proportional to the concentration of uric acid in the sample and is measured using a bichromatic (293,700 nm) endpoint technique.

Reagents

Wells ^a	Form	Ingredient	Concentration ^b	Source
1 – 3, 7	Liquid	Buffer Stabilizers Uricase		
8	Liquid	Stabilizers	8 IU/mL	Bacterial

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal value per test at manufacture.

Precautions: Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact and ingestion.

For *in vitro* diagnostic use

Reagent Preparation: All reagents are liquid and ready to use.

Store at: 2 – 8 °C

Expiration: Refer to carton for expiration date of individual unopened reagent cartridges. Sealed cartridge wells on the instrument are stable for 30 days.

Open Well Stability: 3 days for wells 1 – 3
30 days for well 8

Specimen Collection and Handling: Serum and plasma can be collected using recommended procedures for collection of diagnostic blood specimens by venipuncture.³

Follow the instructions provided with your specimen collection device for use and processing.⁴

Complete clot formation should take place before centrifugation.⁵

Specimens should be free of particulate matter.

Collections of 24-hour urine specimens should be made in containers with NaOH [10 mL of 5% (w/v) NaOH] to prevent urate precipitation.⁶

Separated specimens are stable for 1 day at room temperature, 3 – 5 days at 2 – 8 °C. For longer storage, specimens may be frozen at -20 °C or colder for up to 6 months. Alkaline urine samples are stable at ambient temperature for 3 – 4 days.^{7,8}

Blood collection tubes containing EDTA, lithium heparin, potassium oxalate or sodium fluoride have no measurable effect on this uric acid method.

Corvac® and SST® collection tubes do not affect the URCA method.

Corvac® is a registered trademark of Monoject, Division of Sherwood Medical, St. Louis, MO.
SST® is a registered trademark of Becton-Dickinson, Rutherford, NJ.

Procedure

Materials Provided	
URCA Flex® reagent cartridge, Cat. No. DF77	
Materials Required But Not Provided	
CHEM I Calibrator, Cat. No. DC18A or DC18B	
Quality Control Materials	

Test Steps

Sampling,^c reagent delivery, mixing, processing and printing of results are automatically performed by the Dimension® system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container (if not a primary tube) must contain sufficient quantity to accommodate the sample volume plus the dead volume; precise container filling is not required.

Test Conditions	
Sample Volume	17 µL, (10 µL) ^d
Reagent 1 Volume	132 µL
Reagent 2 Volume	26 µL
Diluent Volume	231 µL
Temperature	37 °C
Wavelength	293 and 700 nm
Type of Measurement	Bichromatic endpoint

d. An alternate sample size of 10 µL can be programmed; refer to the Operator's Guide for the use of alternate sample size.

Calibration

Assay Range	0 – 20.0 mg/dL [0 – 1190 µmol/L] ^e
Calibration Material	CHEM I Calibrator, Cat. No. DC18A or DC18B
Calibration Scheme	3 levels, n = 3
Units	mg/dL [µmol/L] (mg/dL x 59.48) = [µmol/L]
Typical Calibration Levels	0, 12, 23 mg/dL [0, 714, 1368 µmol/L]
Calibration Frequency	Every 3 months for any one lot
A new calibration is required	<ul style="list-style-type: none"> For each new lot of Flex® reagent cartridges After major maintenance or service, if indicated by quality control results As indicated in laboratory quality control procedures When required by government regulations
Assigned Coefficients	Standard sample size = 17 µL C ₀ 0.146 C ₁ -0.070 Alternate sample size = 10 µL C ₀ 0.001 C ₁ -0.123

e. Système International d'Unités [SI Units] are in brackets.

Quality Control

At least once each day of use, analyze two levels of a Quality Control (QC) material with known uric acid concentrations.

Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

Results: The instrument automatically calculates and prints the concentration of uric acid in mg/dL [µmol/L] using the calculation scheme illustrated in your Dimension® Operator's Guide.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Analytical Measurement Range (AMR): 0 – 20.0 mg/dL [0 – 1190 µmol/L]

This is the range of analyte values that can be directly measured on the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

Samples with results in excess of 20.0 mg/dL [1190 µmol/L] should be repeated on dilution.

Manual Dilution: Serum/Plasma: Make appropriate dilution with Reagent grade water to obtain result within the assay range. Enter dilution factor.
Urine: Dilute 1 part urine: 9 parts Reagent grade water. Enter dilution factor of 10.
Serum/Plasma/Urine: Reassay. Resulting readout is corrected for dilution.

Autodilution (AD) (for serum, plasma, urine): Refer to your Dimension® Operator's Guide.
Automated Urine Dilution (AUD) (for Urine): Refer to your Dimension® Operator's Guide.

Limitations of Procedure

The instrument reporting system contains error messages to warn the operator of specific malfunctions. Any report slip containing such error messages should be held for follow-up. Refer to your Dimension® Operator's Guide.

A system malfunction may exist if the following 5-test precision is observed at the standard sample size (17 µL):

Concentration	SD
4.9 mg/dL [291 µmol/L]	>0.2 mg/dL [12 µmol/L]
17.7 mg/dL [1053 µmol/L]	>0.4 mg/dL [24 µmol/L]

Interfering Substances

Xanthine has been reported to decrease the URCA result by 40%.⁹

Formaldehyde (formalin) has been reported to give negative interference with the uricase methods.⁹

Lipemia (Intralipid®) at 600 mg/dL [6.78 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Using the alternate sample size of 10 µL, lipemia (Intralipid®) at 1000 mg/dL [11.3 mmol/L] and above tripped a test report message; therefore the magnitude of the interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

Expected Values

Serum ^d	
Female	2.6 – 6.0 mg/dL [155 – 357 µmol/L]
Male	3.5 – 7.2 mg/dL [208 – 428 µmol/L]
Urine ¹⁰	150 – 990 mg/24hr [0.89 – 5.89 mmol/24hr]

Each laboratory should establish its own reference interval for uric acid as performed on the Dimension® system.

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Specific Performance Characteristics¹

	Precision ^a		
Material	Mean mg/dL [μmol/L]	Standard Deviation (% CV)	
		Within-run	Between-day
Serum			
QCS™ Control Serum Assayed			
Normal	5.1 [303]	0.07 [4] (1.4)	0.07 [4] (1.4)
Abnormal	9.0 [535]	0.11 [6] (1.2)	0.12 [7] (1.3)
Moni-Trol® ^b			
Level 1	3.6 [214]	0.16 [9.5] (4.4)	0.15 [8.9] (4.1)
Level 2	9.7 [577]	0.10 [5.9] (1.1)	0.16 [9.5] (1.6)
Urine			
Fisher Urine Chemistry Control			
Level 1	21.6 [1285]	0.6 [36] (2.8)	0.8 [48] (3.7)
Level 2	27.9 [1660]	0.4 [24] (1.4)	0.7 [42] (2.5)
Urine Pool	37.5 [2230]	0.5 [30] (1.3)	0.8 [48] (2.1)
Urichem Trak ^b			
Level 1	13.6 [809]	0.72 [43] (5.3)	0.88 [52] (6.5)
Level 2 (spiked)	35.4 [2106]	0.80 [48] (2.3)	2.95 [175] (8.4)

f. All specific performance characteristics tests were run after normal recommended equipment quality control checks were performed (refer to your Dimension® Operator's Guide).

g. Specimens at each level were analyzed in triplicate for 20 days. The within-run and between-day standard deviations were calculated by the analysis of variance method.

h. Using reduced sample size (10 μL).

OCS™ is a trademark of Gilford, Irvine, CA 92714.

Moni-Trol® is a registered trademark of Medical Analysis Systems Inc., Camarillo, CA 93012-8058.

Method Comparison

Comparative Method	Regression Statistics			n
	Slope	Intercept mg/dL [μmol/L]	Correlation Coefficient	
Serum URCA method on the aca®				
discrete clinical analyzer ⁱ	0.968	0.16 [9]	0.997	101
Reduced vs. Standard sample size ^j	1.030	-0.03 [-2]	0.999	60 ^k
Urine				
URCA method on the aca® discrete clinical analyzer ⁱ	0.997	-0.69 [-41]	0.962	114

i. Model equation for regression statistics is: Results of Dimension® URCA system = [slope x (comparative method results)] + intercept.

j. Model equation for regression statistics is: Results of Dimension® URCA system using reduced sample size (10 μL) = [slope x results of Dimension® URCA system using standard sample size (17 μL)] + intercept.

k. Range of samples: 0.8 – 14.6 mg/dL [48 – 868 μmol/L].

Specificity

HIL Interference			
The URCA method (using the standard sample size of 17 μL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".			
Substance Tested	Test Concentration [SI Units]	URCA Concentration mg/dL [μmol/L]	Bias ¹ %
Hemoglobin (hemolysate)	1000 mg/dL	6.4 [382]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL [1368 μmol/L]	6.5 [383]	<10
Lipemia (Intralipid®)	200 mg/dL [2.26 mmol/L]	6.3 [377]	<10

The URCA method (using the alternate sample size of 10 μL) was evaluated for interference from hemolysis, icterus and lipemia according to CLSI/NCCLS EP7-P Bias, defined as the difference between the control sample (does not contain interferent) and the test sample (contains the interferent), is shown in the table below. Bias exceeding 10% is considered "interference".

Substance Tested	Test Concentration [SI Units]	URCA Concentration mg/dL [μmol/L]	Bias ¹ %
Hemoglobin (hemolysate)	1000 mg/dL	6.7 [397]	<10
	[0.62 mmol/L] (monomer)		
Bilirubin (unconjugated)	80 mg/dL [1368 μmol/L]	6.7 [397]	<10
Lipemia (Intralipid®)	600 mg/dL [6.78 mmol/L]	6.6 [377]	<10

i. Analyte results should not be corrected based on this bias.

Non-Interfering Substances

The following substances do not interfere with the URCA method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at the URCA concentration of 7.3 mg/dL [434 μmol/L].

Substance	Test Concentration	SI Units
Acetaminophen	0.025 mg/dL	1.66 μmol/L
Amikacin	15 mg/dL	256 μmol/L
Ampicillin	5.3 mg/dL	152 μmol/L
Ascorbic Acid	5 mg/dL	227 μmol/L
Caffeine	6 mg/dL	308 μmol/L
Carbamazepine	3 mg/dL	127 μmol/L
Chloramphenicol	5 mg/dL	155 μmol/L
Chlordiazepoxide	1 mg/dL	33.3 μmol/L
Chlorpromazine	0.2 mg/dL	6.27 μmol/L
Cholesterol	500 mg/dL	12.9 mmol/L
Cimetidine	2 mg/dL	79.2 μmol/L
Creatinine	30 mg/dL	2652 μmol/L
Dextran 40	6000 mg/dL	1500 μmol/L
Diazepam	0.5 mg/dL	17.6 μmol/L
Erythromycin	6 mg/dL	81.6 μmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 μmol/L
Furosemide	6 mg/dL	181 μmol/L
Gentamicin	12 mg/dL	251 μmol/L
Heparin	3 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 μmol/L
Immunoglobulin G	5 g/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 μmol/L
Lithium	2.3 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 μmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 μmol/L
Phenobarbital	10 mg/dL	421 μmol/L
Phenytoin	5 mg/dL	198 μmol/L
Primidone	4 mg/dL	183 μmol/L
Propoxyphene	0.2 mg/dL	4.91 μmol/L
Protein: Albumin	6 g/dL	60 g/L
Protein: Total	12 g/dL	120 g/L
Salicylic Acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 μmol/L
Urea	500 mg/dL	83.3 mmol/L
Valproic Acid	50 mg/dL	3467 μmol/L

Analytical Sensitivity: 0.1 mg/dL [6 μmol/L]


The analytical sensitivity represents the lowest concentration of uric acid that can be distinguished from zero. This sensitivity is defined as the mean value (n = 20) plus two standard deviations of the level 1 (0.0 mg/dL) [0 μmol/L] CHEM I Calibrator.

Symbols Key: See adjacent panel.

Bibliography: See adjacent panel.

Dimension®, Flex® and aca® are trademarks of Siemens Healthcare Diagnostics Inc.

©2009 Siemens Healthcare Diagnostics Inc.
All rights reserved.

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