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The following data represent typical performance for the Dimension<sup>®</sup> clinical chemistry systems and were collected on a Dimension<sup>®</sup> RxL.

Material	Mean U/L [ $\mu$ kat/L]	Precision <sup>7,1</sup>	
		Standard Deviation (%CV)	
		Repeatability	Within-Lab
Multiquel <sup>®</sup> Assayed QC			
Level 1	37 [0.62]	0.6 [0.01] (1.5)	1.6 [0.03] (4.2)
Level 2	157 [2.62]	1.0 [0.02] (0.6)	3.7 [0.06] (2.3)
Level 3	303 [5.06]	3.3 [0.06] (1.1)	7.1 [0.12] (2.4)
Serum Pool 1	81 [1.35]	1.1 [0.02] (1.4)	1.8 [0.03] (2.2)
Serum Pool 2	842 [14.06]	5.7 [0.10] (0.7)	13.3 [0.22] (1.6)

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#### Method Comparison<sup>8</sup> Regression Statistics<sup>9</sup>

Comparative Method	Slope	Intercept U/L [ $\mu$ kat/L]	Correlation Coefficient	n
ADVIA <sup>®</sup> ALPAMP	1.06	-0.4 [-0.01]	0.999	116 <sup>h</sup>

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h. The range of alkaline phosphatase values in the correlation study was 20 – 959 U/L [0.33 – 16.02  $\mu$ kat/L].

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# SIEMENS

REF DF150

**Dimension®** clinical chemistry system

PI-LAB-020/00(01/10/2560)

**Flex® reagent cartridge**

**ALPI**


Issue Date 2012-09-04

### Alkaline Phosphatase

**Intended Use:** The ALPI method is an *in vitro* diagnostic test for the quantitative measurement of alkaline phosphatase in human serum and plasma on the Dimension® clinical chemistry system. Measurements of alkaline phosphatase or its isoenzymes are used in the diagnosis and treatment of liver, bone, parathyroid, and intestinal diseases.

**Summary:** The Dimension® ALPI method is based on the primary reference procedure for the measurement of catalytic activity of alkaline phosphatase at 37 °C as described by the International Federation of Clinical Chemistry (IFCC). The alkaline phosphatase method is based on a procedure published by Bowers and McComb<sup>1</sup> and more recently reviewed by Rej<sup>2</sup>. This method responds to all alkaline phosphatase isoenzymes in human serum.

**Principles of Procedure:** Alkaline phosphatase catalyzes the transphosphorylation of p-nitrophenylphosphate (p-NPP) to p-nitrophenol (p-NP) in the presence of the transphosphorylating buffer, 2-amino-2-methyl-1-propanol (AMP). The reaction is enhanced through the use of magnesium and zinc ions. The change in absorbance at 405 nm due to the formation of p-NP is directly proportional to the ALP activity, since other reactants are present in non-rate limiting quantities and is measured using a bichromatic (405, 510 nm) rate technique.

Reagents Wells <sup>a</sup>	Form	Ingredient	Concentration <sup>b</sup>
1 – 6 (Reagent 1)	Liquid	2-Amino-2-Methyl-1-Propanol (AMP) Magnesium Acetate Zinc Sulfate HEDTA	3.0 M 8.0 mmol/L 4.0 mmol/L 8.0 mmol/L
7 – 8 (Reagent 2)	Liquid	p-NPP Buffer	101.6 mmol/L

a. Wells are numbered consecutively from the wide end of the cartridge.

b. Nominal final value per well in a cartridge.

### Risk and Safety:

Irritant. Contains 2-amino-2-methyl-1-propanol.  
R36/38: Irritating to eyes and skin.

R52/53: Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Safety data sheets (MSDS/SDS) available on [www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

### Precautions:

Used cuvettes contain human body fluids; handle with appropriate care to avoid skin contact or ingestion.

For *in vitro* diagnostic use.

**Reagent Preparation:** All reagents are liquid and ready to use.

**Store at:** 2 – 8 °C

Protect from light after opening.

**Expiration:** Refer to carton for expiration date of individual unopened reagent cartridges. Sealed wells on the instrument are stable for 30 days.

**Open Well Stability:** 2 days for wells 1 – 6  
4 days for wells 7 – 8

### Specimen Collection and Handling:

Recommended specimen types: serum and plasma (lithium heparin).

Normal procedures for collecting and storing serum and plasma may be used for samples to be analyzed by this method. In the preparation of serum or plasma samples, avoid prolonged contact with separated red cells.<sup>3</sup>

Follow the instructions provided with your specimen collection device for use and processing.<sup>4</sup>

For Serum: Complete clot formation should take place before centrifugation. Serum or plasma should be physically separated from cells as soon as possible with a maximum limit of two hours from the time of collection.<sup>5</sup>

Specimens are stable for 8 hours at room temperature, 7 days at 2 – 8 °C and 6 months when frozen at -20 °C or colder.<sup>6</sup> Avoid repeated freezing and thawing. Thawed specimens which are turbid must be clarified by centrifugation prior to testing.

The purpose of specimen handling and storage information is to provide guidance to users; however, users may validate their own procedures for handling and storing patient samples.

### Procedure

#### Materials Provided

ALPI Flex® reagent cartridge, Cat. No. DF150

#### Materials Required But Not Provided

ALPI CAL, Cat. No. DC150

Enzyme Diluent, Cat. No. 790035901

Quality Control Materials

#### Test Steps

Sampling,<sup>7</sup> reagent delivery, mixing, and processing are automatically performed by the Dimension® clinical chemistry system. For details of this processing, refer to your Dimension® Operator's Guide.

c. The sample container must contain sufficient quantity to accommodate the sample volume plus dead volume. Precise container filling is not required.

#### Test Conditions

Sample Volume (delivered to the cuvette)	7 µL
Reagent 1 Volume	90 µL
Reagent 2 Volume	57 µL
Diluent Volume	206 µL
Temperature	37 °C
Reaction Time	7.2 minutes
Wavelength	405 and 510 nm
Type of Measurement	Bichromatic rate

#### Calibration

Assay Range	10 – 1000 U/L [0.17 – 16.70 µkat/L] <sup>d</sup>
Calibration Material	Dimension® ALPI CAL, Cat. No. DC150
Calibration Scheme	Three levels (n = 3)
Units	U/L [µkat/L] (U/L x 0.0167) = [µkat/L]
Typical Calibration Levels	Level 1: 0 U/L [0.00 µkat/L] Level 2: 500 U/L [8.35 µkat/L] Level 3: 1000 U/L [16.70 µkat/L]
Calibration Frequency	Every 90 days for any one lot
A new calibration is required:	<ul style="list-style-type: none"> <li>For each new lot of Flex® reagent cartridges</li> <li>After major maintenance or service, if indicated by quality control results</li> <li>As indicated in laboratory quality control procedures</li> <li>When required by government regulations</li> </ul>
Assigned Coefficients	C <sub>0</sub> 0.000 C <sub>1</sub> 1.000

d. Système International d'Unités [SI Units] are in brackets.

#### Quality Control

Follow government regulations or accreditation requirements for quality control frequency. At least once each day of use, analyze two levels of a Quality Control (QC) material with known alkaline phosphatase activity. Follow your laboratory internal QC procedures if the results obtained are outside acceptable limits.

**Results:** The instrument calculates the activity of alkaline phosphatase in U/L [µkat/L] using the calculation scheme described in your Dimension® Operator's Guide.

**Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.**

#### Analytical Measurement Range (AMR): 10 – 1000 U/L [0.17 – 16.70 µkat/L]

This is the range of analyte values that can be measured directly from the specimen without any dilution or pretreatment that is not part of the usual analytical process and is equivalent to the assay range.

• Samples with results in excess of 1000 U/L [16.70 µkat/L] are reported as "Above Assay Range" and should be repeated on dilution.

**Autodilution (AD):** The autodilute sample volume is 3 µL (dilution factor = 2.3) for serum and plasma. Refer to your Dimension® Operator's Guide.

**Manual Dilution:** Results in excess of 1000 U/L [16.70 µkat/L] should be repeated after diluting the sample with Enzyme Diluent (Cat. No. 790035901) or equivalent to produce a sample within the assay range.<sup>8</sup> Enter dilution factor. Reassay. Resulting readout is corrected for dilution.

• Samples with results less than 10 U/L [0.17 µkat/L] will be reported as "less than 10 U/L [0.17 µkat/L]."

e. Specimens should not be diluted to a final activity less than 200 U/L [3.34 µkat/L]

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#### Limitations of Procedure

The instrument reporting system contains flags and comments to provide the user with information regarding instrument processing errors, instrument status information and potential errors in alkaline phosphatase results. Refer to your Dimension® Operator's Guide for the meaning of report flags and comments. Any report containing flags and/or comments should be addressed according to your laboratory's procedure manual.

A system malfunction may exist if the following 5-test precision is observed:

Activity	SD
150 U/L [2.51 µkat/L]	> 7 U/L [0.12 µkat/L]
500 U/L [8.35 µkat/L]	> 16 U/L [0.27 µkat/L]

#### Interfering substances

Triglycerides above 1500 mg/dL [16.95 mmol/L] tripped a test report message; interference could not be determined.

Lipemia (Intralipid®) at 600 mg/dL [6 g/L] and above tripped a test report message; interference could not be determined.

Intralipid® is a registered trademark of Fresenius Kabi AG, Bad Homburg, Germany.

**Expected Values: 46 – 116 U/L [0.77 – 1.90 µkat/L]**

The reference interval was calculated non-parametrically and represents the central 95% of results determined from a population of healthy adults (n = 132).

Each laboratory should establish its own expected values for alkaline phosphatase as performed on the Dimension® clinical chemistry system.

#### Specific Performance Characteristics

The following data represent typical performance for the Dimension® clinical chemistry systems and were collected on a Dimension® RxL.

Material	Mean U/L [µkat/L]	Precision <sup>2,1</sup>	
		Standard Deviation (%CV)	Repeatability
			Within-Lab
Multiquel® Assayed QC			
Level 1	37 [0.62]	0.6 [0.01] (1.5)	1.6 [0.03] (4.2)
Level 2	157 [2.62]	1.0 [0.02] (0.6)	3.7 [0.06] (2.3)
Level 3	303 [5.06]	3.3 [0.06] (1.1)	7.1 [0.12] (2.4)
Serum Pool 1	81 [1.35]	1.1 [0.02] (1.4)	1.8 [0.03] (2.2)
Serum Pool 2	842 [14.06]	5.7 [0.10] (0.7)	13.3 [0.22] (1.6)

f. CLSI/NCCLS EP5-A2 was used. During each day of testing, two separate runs, with two test samples, for each test material, were analyzed for 20 days.

Multiquel® Assayed Control is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

Comparative Method	Slope	Method Comparison <sup>3</sup> Regression Statistics <sup>4</sup>		n
		Intercept U/L [µkat/L]	Correlation Coefficient	
ADVIA® ALPAMP	1.06	-0.4 [-0.01]	0.999	116 <sup>5</sup>

g. CLSI/NCCLS EP9-A2 was used. The method used to fit the linear regression line was ordinary least squares.

h. The range of alkaline phosphatase values in the correlation study was 20 – 959 U/L [0.33 – 16.02 µkat/L].

#### Serum and Plasma Equivalency

No clinically significant difference was observed between serum (x) and plasma (y) as shown in the following ordinary least squares regression.

Sample Comparison	Slope	Intercept U/L [µkat/L]	Correlation Coefficient	n
Lithium heparin plasma vs. Serum	1.02	-5 [-0.08]	0.999	50 <sup>6</sup>

i. The range of alkaline phosphatase values in the correlation study was 39 – 884 U/L [0.65 – 14.76 µkat/L].

#### Specificity

##### Hemolysis, Icterus, Lipemia (HIL) Interference

The ALPI method was evaluated for interference according to CLSI/NCCLS EP7-A2.<sup>8</sup> Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Substance Tested	Substance Concentration	Alkaline Phosphatase U/L [µkat/L]	Bias* %
Hemoglobin (hemolysate)	1000 mg/dL [0.62 mmol/L]	297 [4.96]	<10
		823 [13.74]	<10
		300 [5.01]	<10
Bilirubin (unconjugated)	80 mg/dL [1026 µmol/L]	834 [13.92]	<10
		307 [5.12]	<10
		856 [14.29]	<10
Bilirubin (conjugated)	80 mg/dL [1026 µmol/L]	308 [5.14]	<10
		878 [14.66]	<10
		308 [5.14]	<10
Lipemia (Intralipid®)	500 mg/dL [5 g/L]	878 [14.66]	<10
		308 [5.14]	<10
		878 [14.66]	<10

\* Analyte results should not be corrected based on this bias.

j. Lipemia (Intralipid®) at 600 mg/dL and above tripped a test report message; interference could not be determined.

#### Non-Interfering Substances

The following substances do not interfere with the ALPI method when present in serum and plasma at the concentrations indicated. Inaccuracies (biases) due to these substances are less than 10% at alkaline phosphatase activity of 299 U/L [4.99 µkat/L] and 842 U/L [14.06 µkat/L].

Substance	Test Concentration	SI Units
Acetaminophen	20 mg/dL	1324 µmol/L
Amikacin	8 mg/dL	137 µmol/L
Ampicillin	5.3 mg/dL	152 µmol/L
Ascorbic acid	6 mg/dL	342 µmol/L
Caffeine	6 mg/dL	308 µmol/L
Carbamazepine	3 mg/dL	127 µmol/L
Chloramphenicol	5 mg/dL	155 µmol/L
Chlordiazepoxide	1 mg/dL	33.3 µmol/L
Chlorpromazine	0.2 mg/dL	6.27 µmol/L
Cholesterol	503 mg/dL	13 mmol/L
Cimetidine	2 mg/dL	79.2 µmol/L
Creatinine	30 mg/dL	2.65 mmol/L
Dextran 40	6000 mg/dL	1500 µmol/L
Diazepam	0.5 mg/dL	18 µmol/L
Digoxin	6.1 ng/mL	7.8 nmol/L
Erythromycin	6 mg/dL	81.6 µmol/L
Ethanol	400 mg/dL	86.8 mmol/L
Ethosuximide	25 mg/dL	1770 µmol/L
Furosemide	6 mg/dL	181 µmol/L
Gentamicin	1 mg/dL	21 µmol/L
Heparin	3.0 U/mL	3000 U/L
Ibuprofen	50 mg/dL	2425 µmol/L
Immunoglobulin G (IgG)	5000 mg/dL	50 g/L
Lidocaine	1.2 mg/dL	51.2 µmol/L
Lithium	2.2 mg/dL	3.2 mmol/L
Nicotine	0.1 mg/dL	6.2 µmol/L
Penicillin G	25 U/mL	25000 U/L
Pentobarbital	8 mg/dL	354 µmol/L
Phenobarbital	10 mg/dL	431 µmol/L
Phenytoin	5 mg/dL	198 µmol/L
Primidone	4 mg/dL	183 µmol/L
Propoxyphene	0.16 mg/dL	4.91 µmol/L
Protein: Total	12000 mg/dL	120 g/L
Protein: Albumin	6000 mg/dL	60 g/L
Salicylic acid	60 mg/dL	4.34 mmol/L
Theophylline	4 mg/dL	222 µmol/L
Triglycerides	1500 mg/dL	16.95 mmol/L
Urea	500 mg/dL	83 mmol/L
Uric acid	20 mg/dL	1.2 mmol/L
Vancomycin	10 mg/dL	69 µmol/L
Valproic acid	50 mg/dL	3467 µmol/L

#### Limit of Detection and Limit of Blank<sup>10,8</sup>

The Limit of Detection (LoD) for ALPI is 8 U/L [0.13 µkat/L] determined consistent with CLSI guideline EP17-A<sup>10</sup> and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 120 determinations, with 5 blank (12 reps per sample) and 5 low level samples (12 reps per sample). The Limit of Blank (LoB) is 3 U/L [0.05 µkat/L].

k. LoD is the lowest concentration of analyte that can be detected reliably. LoB is the highest concentration that is likely to be observed for a blank sample.

#### Limit of Quantitation<sup>10,9</sup>

The Limit of Quantitation (LoQ) for ALPI is 10 U/L [0.17 µkat/L] determined consistent with CLSI Guideline EP17-A and allowable total error of 40%.

l. LoQ is the lowest amount of alkaline phosphatase that can be quantitatively determined within a defined total error.

**Symbols Key:** See adjacent panel.

**Bibliography:** See adjacent panel.

Dimension®, ADVIA® and Flex® are trademarks of Siemens Healthcare Diagnostics.

© 2012 Siemens Healthcare Diagnostics

All rights reserved.

[illegible]